

AD-A119 041

AIR FORCE INST OF TECH WRIGHT-PATTERSON AFB OH
COMPUTER AUTOMATION OF THE THERMAL PULSE TECHNIQUE FOR LOCAL BL--ETC(U)
1982 K L BAUM
AFIT/CI/NR/82-45T

F/G 6/16

NL

UNCLASSIFIED

[or]

AD A
13041



END
DATE
FILMED
10-82
DTIC

UNCLASS

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

REPORT DOCUMENTATION PAGE		READ INSTRUCTIONS BEFORE COMPLETING FORM
1. REPORT NUMBER AFIT/CI/NR/82-45T	2. GOVT ACCESSION NO.	3. RECIPIENT'S CATALOG NUMBER
4. TITLE (and Subtitle) Computer Automation of the Thermal Pulse Technique for Local Blood Flow Measurements		5. TYPE OF REPORT & PERIOD COVERED THESIS/DISSERTATION/
7. AUTHOR(s) Kurt Lewis Baum		6. PERFORMING ORG. REPORT NUMBER
9. PERFORMING ORGANIZATION NAME AND ADDRESS AFIT STUDENT AT: University of Illinois		8. CONTRACT OR GRANT NUMBER(s)
11. CONTROLLING OFFICE NAME AND ADDRESS AFIT/NR WPAFB OH 45433		10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office)		12. REPORT DATE 1982
		13. NUMBER OF PAGES 79
		15. SECURITY CLASS. (of this report) UNCLASS
		15a. DECLASSIFICATION/DOWNGRADING SCHEDULE
16. DISTRIBUTION STATEMENT (of this Report) APPROVED FOR PUBLIC RELEASE; DISTRIBUTION UNLIMITED		
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)		
18. SUPPLEMENTARY NOTES APPROVED FOR PUBLIC RELEASE: IAW AFR 190-17 30 AUG 1982 LYNN E. WOLAVER Dean for Research and Professional Development AFIT, Wright-Patterson AFB OH		
19. KEY WORDS (Continue on reverse side if necessary and identify by block number)		
20. ABSTRACT (Continue on reverse side if necessary and identify by block number) ATTACHED		

AD A119041

FILE COPY

DD FORM 1473
1 JAN 73

EDITION OF 1 NOV 65 IS OBSOLETE

UNCLASS

82 09 07 418

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

COMPUTER AUTOMATION OF THE THERMAL PULSE TECHNIQUE FOR
LOCAL BLOOD FLOW MEASUREMENTS. Kurt Lewis Baum, USAF, 1982,
79 pages, Master of Science in Electrical Engineering, University
of Illinois.

Tissue blood perfusion is a fundamental measurement in
physiology that affects the entire spectrum of medical practice
and research. A new and innovative method is under development
by Dr. Kenneth R. Holmes and Dr. Michael M. Chen at the
University of Illinois. Their thermal pulse-decay method
utilizes a small thermistor to pulse heat the tissue under
study. The thermistor is then used to record tissue temperature
as the heat dissipates due to thermal conductivity and blood
perfusion. From this cooling data, local blood perfusion can
be calculated by various computer routines.

The process of initiating and controlling the experiment,
acquiring and storing the data, and calculating perfusion
parameters has been computer automated. The system is based
on a Digital Equipment Corporation LSI 11 minicomputer. The
software package developed for the system is user oriented.
It can control up to six probes at once, performing both
heating and measurement tasks. The user is free to choose
the duration of the heat pulses, as well as the sampling rate
and sampling duration after the heat pulse. The program
automatically generates a data file for each active probe.
The files can be recalled for display on a graphics terminal
or for calculating perfusion parameters. In addition, an
automatic mode is available which repetitiously performs the
experiment with no user interaction.

The computer automation of the thermal pulse technique
for local blood flow measurements will allow further development
of this promising new measurement tool.

Accession For	
NTIS GRA&I	<input checked="checked" type="checkbox"/>
DTIC TAB	<input type="checkbox"/>
Unannounced	<input type="checkbox"/>
Justification	
By	
Distribution/	
Availability Codes	
Dist	Avail and/or Special

11



AFIT/CI/NR/82-45T

AFIT/NR
Wright-Patterson AFB OH 45433

AUTHOR: Kurt Lewis Baum

1. Did this research contribute to a current Air Force project?

() b. NO

() a. YES

() b. NO

() a. MAN-YEARS

() b. \$

() a. HIGHLY
SIGNIFICANT

() b. SIGNIFICANT

() c. SLIGHTLY

() d. OF NO
SIGNIFICANCE

5. AFIT welcomes any further comments you may have on the above questions, or any additional details concerning the current application, future potential, or other value of this research. Please use the bottom part of this questionnaire for your statement(s).

POSITION

LOCATION

STATEMENT(s):

COMPUTER AUTOMATION OF THE THERMAL PULSE
TECHNIQUE FOR LOCAL BLOOD FLOW MEASUREMENTS

BY

KURT LEWIS BAUM

B.S., United States Air Force Academy, 1981

THESIS

Submitted in partial fulfillment of the requirements
for the degree of Master of Science in Electrical Engineering
in the Graduate College of the
University of Illinois at Urbana-Champaign, 1982

Urbana, Illinois

COMPUTER AUTOMATION OF THE THERMAL PULSE TECHNIQUE FOR LOCAL BLOOD FLOW MEASUREMENTS. Kurt Lewis Baum, USAF, 1982, 79 pages, Master of Science in Electrical Engineering, University of Illinois.

Tissue blood perfusion is a fundamental measurement in physiology that affects the entire spectrum of medical practice and research. A new and innovative method is under development by Dr. Kenneth R. Holmes and Dr. Michael M. Chen at the University of Illinois. Their thermal pulse-decay method utilizes a small thermistor to pulse heat the tissue under study. The thermistor is then used to record tissue temperature as the heat dissipates due to thermal conductivity and blood perfusion. From this cooling data, local blood perfusion can be calculated by various computer routines.

The process of initiating and controlling the experiment, acquiring and storing the data, and calculating perfusion parameters has been computer automated. The system is based on a Digital Equipment Corporation LSI 11 minicomputer. The software package developed for the system is user oriented. It can control up to six probes at once, performing both heating and measurement tasks. The user is free to choose the duration of the heat pulses, as well as the sampling rate and sampling duration after the heat pulse. The program automatically generates a data file for each active probe. The files can be recalled for display on a graphics terminal or for calculating perfusion parameters. In addition, an automatic mode is available which repetitiously performs the experiment with no user interaction.

The computer automation of the thermal pulse technique for local blood flow measurements will allow further development of this promising new measurement tool.

ACKNOWLEDGEMENTS

I wish to express my sincere appreciation to my advisor, Dr. Richard L. Magin, for providing me the opportunity to work on this project. I am also grateful for his support and guidance throughout the past year.

A very special thanks is extended to Dr. Kenneth R. Holmes, without whom this project would not have been possible. His research goals and desires were the determining factors in the direction of this project.

TABLE OF CONTENTS

<u>CHAPTER</u>	<u>Page</u>
I. INTRODUCTION.	1
II. THERMAL PULSE-DECAY METHOD.	6
III. HARDWARE.	11
IV. SOFTWARE.	17
V. RESULTS	27
VI. CONCLUSIONS AND RECOMMENDATIONS	36
APPENDIX A. SYSTEM HARDWARE CONNECTIONS.	37
APPENDIX B. DATA FILE STRUCTURES	38
APPENDIX C. OPERATING PROGRAM.	40

I. INTRODUCTION

Background

Tissue blood perfusion is a fundamental measurement in physiology that affects the entire spectrum of medical practice and research [1]. Many different techniques have been developed to measure tissue blood flow. Two common methods involve indicator dilution and radio-labeled microspheres [2]. These measurements are complicated and cannot be repeated at frequent intervals.

A different category of blood perfusion measurements involves the use of thermal techniques. These methods have the potential to overcome the limitations of the indicator dilution and microsphere techniques. However, most thermal methods have the drawback of comparing the heat dissipation against a value of thermal conductivity for nonperfused tissue. This requires either a suspension of the blood flow to the tissue or the use of tabulated values for the type of tissue being examined. Suspension of blood flow is traumatic to tissue and might mean sacrificing the animal, while using tabulated values can lead to inaccuracies due to assumptions made in creating the table. In addition, the theoretical basis for some of these methods is open to question, since the volume of tissue being sampled is not much larger than the probe itself, and probe size and shape have been shown to alter results [1].

The thermal pulse-decay method being developed at the University of Illinois by Dr. Kenneth R. Holmes and Dr. Michael M. Chen is a thermal method that overcomes all of the problems described above. This method consists of inserting a small thermistor into the tissue of interest. A known quantity of heat is deposited in the tissue when current flows through the thermistor. The thermistor is then used to measure post-pulse tissue cooling, from which local perfusion and thermal conductivity can be determined [1].

The basis of these calculations arises from the heat dissipation mechanisms working in the tissue. With no blood perfusion, the primary decay in temperature is due to thermal diffusion. This temperature decay can be mathematically modeled as a decaying power series [1]. Blood perfusion is a thermal transfer mechanism that can be modeled as an exponential decay of temperature. In vivo tissue cooling typically includes both of these mechanisms and can be modeled as a product of the exponential and power series. From the shape of this curve, the thermal conductivity and the local blood perfusion can be determined [1].

The advantages of this method are: it provides an absolute measurement of the volumetric perfusion rate (ml blood/ml tissue sec) without requiring calibrations or stop-flow measurements; the sampling volume is considerably larger than the volume of tissue traumatized by the microprobe; the probe shape and size do not affect the results; the electronics and

calculations are extremely simple; the increase in tissue temperature is usually only 0.5 C; the small diameter of the probe causes minimal trauma in the tissue under examination [3,4].

Some drawbacks of this method are that it does require insertion of a probe, it yields point information, and it is an indirect measurement process.

Problem and Scope

Small thermistor beads are fabricated into needle-like probes to aid in insertion into the tissue and minimization of trauma in the tissue. The thermistor is incorporated into a bridge circuit that performs both heating and measurement roles. The bridge output is used to generate cooling curves on a strip chart recorder. Computer programs have been written to analyse these cooling curves for thermal conductivity and perfusion rate.

The problem with this procedure is that the data must be read point by point from the the cooling curves and typed into the computer. This process is both time consuming and inaccurate. Additionally, since the data is analyzed after the experiment has been completed, the operator has no opportunity to modify the experimental parameters.

The solution to these problems is to automate the control of the experiment, the acquisition and storage of the data, and calculations performed on the data. This report describes a

computer based hardware and software system that will not only solve the problems of entering cooling curves into the computer and freeing the operator from continuous supervision of the system, but should also give the operator the advantage of having real-time blood flow data to tailor each succeeding measurement. Figure 1 is a block diagram of the proposed system.

Presentation

Chapter II describes the characteristics of the thermistor bridge circuit, the governing equations, and the desired specifications for the proper operation of the control system.

Chapter III describes the hardware used by the system while Chapter IV documents the system software.

Chapter V presents the results of experimental testing of the entire system.

Chapter VI contains the conclusions about the present system and recommendations for future systems.

Software flowcharts and a complete program listing are contained in the appendices.




Figure 1. System Block Diagram

II. THERMAL PULSE-DECAY METHOD

Description

The bridge circuit containing the thermistor probe is a very simple balancing circuit with one modification that allows current to pass through the probe during heating. This circuit is shown in Figure 2. The balance resistor, R_b , is used to adjust the bridge output. A five Volt signal impressed at V_t can drive the transistor switch for the heater relay. The duration of this signal determines the length of the heating pulse and thus the energy delivered to the tissue.

Simple circuit analysis of the bridge yields the following relationship between the output voltage, E , and the probe resistance, R_p .


$$R_p = \left(\frac{R_s V}{R_s V / (R_b + R_s) + E} \right) - R_s \quad (1)$$

Then, once the resistance of the thermistor is known, the temperature, T , can be found as:

$$T = A - B \ln(R_p) \quad (2)$$

where A and B are calibrated constants of the particular thermistor probe in use [1].

The temperature rise due to tissue heating causes a corresponding increase in the output voltage. Thus, even though a nonlinear relationship exists between voltage and



- $V_h = 4.23$ Volts
- $V_b = 5.45$ Volts
- $V_r = 12$ Volts
- $R_p =$ Thermistor
- $R_b =$ Bridge Balance Potentiometer
(1000 to 1500 Ohms)
- $R_s = 22.05$ kOhms
- $R_t = 3.3$ kOhms

Figure 2. Thermistor Bridge Circuit

temperature, a cooling curve obtained by plotting the output voltage of the bridge looks similar to a cooling curve that plots actual temperature.

For a cooling curve expressed in Volts, the typical experiment will have about a seven millivolt range. In addition to the actual cooling part of the curve, important information is contained in the temperature values immediately prior to heating. This pre-sample period can be used to project a baseline for use in normalizing the cooling part of the curve to temperature drifts. A sample cooling curve is shown in Figure 3.

System Specifications

The following is a list of desired specifications for the automated control and data acquisition system.

1. Deliver a square wave five Volt pulse to each of six bridge circuits to control the heating cycle.
2. Vary the duration of each pulse length independently over a range of zero to twenty seconds in one-tenth second intervals.
3. The duration time of each heat pulse must be stored for computation purposes.
4. Monitor and record voltage output of six thermistor bridges with a resolution of three millivolts and a range of minus five to five volts.
5. Allow for sampling rates up to sixty Hertz in each of six channels.




Figure 3. Typical Cooling Curve

6. Channel selection may be varied without loss of information. Channels need not be operated in sequential order.
7. Stored data may be recalled for display on a video terminal or printed on hard copy.
8. System must include a file management system for stored data files.
9. Calculations may be performed while system is not being used to measure perfusion.
10. Calculate and display steady state body temperature as calculated from thermistor data. Calibrate system for each thermistor.
11. If input voltage exceeds an operator selected value, system will issue an audible alarm. Alarm may be disabled.
12. Repeat measurements (delivery of pulses) automatically.

III. HARDWARE

The system hardware is composed of three main units. These are the thermistor bridge, the signal conditioner, and the computer. Figure 4 is a block diagram of the entire computer system, with the computer separated into its major components.

Bridge Circuit

The bridge circuit is described in Chapter II.

Signal Conditioning

Preliminary tests of the bridge circuit showed considerable noise to be present, which appeared in two basic forms. The most prevalent noise was sixty cycle power line interference, which had a magnitude up to four millivolts, almost as large as the desired signal. The second kind of noise was a very high frequency noise whose source could not be identified.

In addition to the noise, the signal strength was on the order of millivolts, too weak to be used as an input signal to the analog to digital converter, which measures in volts. The signal conditioning block, shown in Figure 5, is used to correct for these two problems. First, the signal was amplified by a differential laboratory amplifier with a gain of 1000. The input signal for a typical experiment has a range of about seven millivolts, with both positive and




Figure 4. System Block Diagram




Figure 5. Signal Conditioning Schematic Diagram

negative components. The Analog to Digital Converter (ADC) has an input range of minus five to plus five volts. If the bridge is balanced at the beginning of an experiment for an output of about three millivolts, the gain of 1000 will provide a signal at the input of the ADC that will range from minus three to plus four volts. This is entirely within the range of the ADC and still allows for some temperature drift in between heating pulses.

A two stage filter was built to eliminate the noise problem. The first stage consists of a second-order band reject filter, the output of which provided 38 decibels of suppression centered at 60.4 Hertz with a quality factor of five. This filter effectively attenuates the sixty cycle interference. The second stage of the filter is a simple resistor-capacitor low pass filter with a 3 dB point of 10.6 Hertz. This filter eliminates the high frequency noise and further reduces the sixty cycle noise. The frequency content of the input signal is much less than 10 Hertz, so the filter does not attenuate the thermal pulse data signal.

Computer

This section will describe the capabilities and functions of each block of the computer. The computer is a Digital Equipment Corporation (DEC) LSI 11 minicomputer. Hardware specifications for the LSI 11 can be found in References [5] and [6].

The Parallel Input-Output unit is a DEC DRV11. It consists of 16 separate input and output lines, along with appropriate control lines. In this application the unit is used to control the heat pulses and light appropriate indicator lights on the bridge unit. In future versions, it can be used to monitor switch positions on the bridge unit. The pin assignments for this board are contained in Appendix A. More information on the operation of this board is contained in Reference [6].

The Analog to Digital converter is an ADAC 1030 which has eight differential inputs each with a range of minus five to plus five volts and a programmable gain of 1, 2, 5, or 10. In future versions of the bridge, this gain can be used as part of the 1000 gain of the signal conditioner. With the full gain of 1000, the ADC's 12 bit resolution is able to look at signals from minus five millivolts to five millivolts in 2.5 microvolt increments. This meets the original specification for resolution. Timing for the sampling of the ADC is accomplished with the real time clock of the LSI 11 computer. This allows sampling rates as high as sixty hertz for each of six channels. Additional information on the operation of the ADC can be found in Reference [7].

The Central Processing Unit is a DEC LSI 11/23 CPU with Memory Management Unit (MMU), model KDF11-AA, along with the floating point hardware option, model KEF11-AA. The MMU will allow full and efficient use of the Chrislin memory board.

The speed and power of the LSI 11/23 CPU with floating point should allow close to real time calculation of the desired parameters after the data has been collected.

The Chrislin CI-1123 memory provides a full 256 kilobytes of random access memory. In addition to containing the operating program, it allows for storage of quite large data arrays during sampling that can be transferred later to disk.

The Data Systems Design 470 disk drive provides two megabytes of on line storage. Drive one will be used as a system disk containing the source program and monitor, while drive two will be used to interchange data disks.

The serial interface is a DEC DLV11-J, which has four independent communication channels. One of these is used for the system terminal. The others could be used for a modem, a remote terminal, or a line printer.

The system terminal is a Lear Siglar ADM5 with a 512 retrographics board. This provides graphics capability that is software compatible to the Techtronics 4010 terminal. The main use of the graphics will be to display various cooling curves for visual inspection. Additional information is contained in References [8] and [9].

The GP100 graphics plotter is connected directly to the ADM5 terminal. It provides hardcopy for any of the graphics sent to the terminal or functions as a line printer for listing programs or data sets. It is further described in Reference [10].

IV. SOFTWARE

The system software is written in FORTRAN and resides on the system disk. Both the user and the operating program have access to an extensive library of programs written for the DEC LSI 11 computer system. These include the RT11 monitor, a text editor, a disk handler, linkers, compilers, various input-output routines, and a system library of FORTRAN callable subroutines. A full description of these programs is contained in Reference [11].

The FORTRAN operating program contained on the system disk allows the user to interact with the system to configure the experiment and data collection in an easy and flexible manner. In addition, it allows calculation routines to be performed on the data and several modes of automatic operation. The user also has the ability to list or graph any data file.

This chapter contains a complete description of the operating program. Figure 6 is a simplified flowchart of the main program. Flowcharts of the subroutines and a complete listing of the FORTRAN code are contained in Appendix C.

Main Program

The program begins with an initialization routine that sets each variable to a valid starting value. As the program is used, the initial values can be modified to reflect more accurately the desired initial conditions for the experiments.




Figure 6. Simplified Flowchart of the Main Operating Program

The terminal is then filled with the current values of the program variables and a list of the commands available to the user. An example of this is shown in Figure 7. The commands are single letters that correspond to the function or parameter to its right. Each second, the current time is displayed, all six bridge outputs are sampled and displayed in both millivolts and degrees Celsius, and the time until the next automatic operation is updated and displayed. While in this mode, the software continuously checks for user input or the conditions necessary to cause an automatic operation.

Each user input causes the main program to call an appropriate subroutine which actually performs the command. A summary of each of the available commands is shown in Table I. The following is a description of each command along with its use and restrictions.

Command 'O'

This command allows the user to choose which channels will be active when either SAMPLE or CALCULATE is performed. Each channel can be independently turned on or off with no restrictions as to which channels must be used for a given number of probes. The software for command 'O' is contained in subroutine UPDATE.

Command 'F'

This command allows the filename of any channel to be changed. Each filename must be constructed of three

```

TIME 13: 3: 16 DATE 5/16/82 DATA VERSION ONE

  CHANNEL          1          2          3          4          5          6
O  SAMPLE STATUS    ON        ON
  CALCULATE STATUS  ON
F  FILENAME         AAA000    AAB000    AAC000    AAD000    AAE000    AAF000
P  PROBE            24        25        0         0         0         0
  CURRENT MILLIVOLTS 0.003    0.100    0.944    1.529    2.200    5.250
  CURRENT TEMPERATURE 34.432    35.585    0.000    0.000    0.000    0.000
  LAST TEMP         34.432    35.022    0.000    0.000    0.000    0.000
H  HEAT (SEC:TICKS)  2: 0     2: 0     2: 0     2: 0     2: 0     2: 0
  DURATION (SEC:TIC) 3: 0     60      60      100     100     1
R  PRESAMPLE        3: 0     60      60      900     900     1
  SAMPLE            15: 0    60      60      900     900     1

L  UPPER ALARM LIMIT (MV) 6.25    A  AUTOMATIC RUN
  LOWER ALARM LIMIT (MV) -6.00    B  BREAK FROM AUTO
E  EXPERIMENT INTERVAL (MIN:SEC) 0:45    C  CALCULATE
  EXPERIMENT REPETITIONS 5          D  LOOK AT DATA FILE
  NEXT EXPERIMENT (MIN:SEC) 0: 0     S  SAMPLE
                                Z  QUIT

PLEASE ENTER YOUR COMMAND

```

Figure 7. Video Terminal Presentation of Program Status and Commands

Table 1. Summary of Commands

<u>Command</u>	<u>Description</u>
O	Change Channel Activity Status
F	Change Data Filename
P	Change Probe, Bridge Balance, or Descriptive Text
H	Change Heat Pulse Duration
R	Change Duration and Frequency of Sampling
L	Change Alarm Limits
E	Change Experiment Repetition Parameters
A	Enter Automatic Operation
B	Break From Automatic Operation
C	Perform Calculations
D	Display Data Files
S	Perform Heat Pulse and Sampling
Z	Exit Program

alphanumeric characters followed by three numeric characters. This allows the user to describe the experiment using three letters and then determine the experiment repetition using the three numerals remaining. This filename is used by both SAMPLE and CALCULATE. Subroutine UPDATE performs the command 'F'.

Command 'P'

This command allows the user to change the probe number, probe calibration data, bridge balance conditions, and descriptive text of any channel. When a probe is changed, the software reads in the calibration data for the new probe from a disk file named PROBE.DAT. This file is maintained by the auxiliary program PROBEC.FOR which is described later. The calibration data is part of the data file stored when a sample is performed. It is needed to calculate absolute temperature.

Changing the bridge balance conditions requires the user to enter the new resistance of the potentiometer. This value is also needed for absolute temperature calculation, and is part of the data file stored by SAMPLE.

The descriptive text can be used to record probe placement, experiment objective or other pertinent information. It is also part of the data file stored by SAMPLE.

Subroutine UPDATE performs the command 'P'.

Command 'H'

This command allows the user to modify the duration of the heat pulse that is applied to the probe. Each channel can be varied independently in increments of one-sixtieth of a second. A duration of zero is allowed and simply means that that channel will not be pulsed during sampling. This allows probes to be used to determine heat patterns delivered by other probes. Subroutine UPDATE performs the command 'H'. The heat pulse duration is part of the data file stored by SAMPLE.

Command 'R'

This command allows the user to define the sampling rate and duration for an experiment. The samples taken before the heat pulse, the presample period, can be configured independently of the samples taken after the heat pulse. The user inputs the desired time duration and frequency of sampling. The program calculates and displays the actual number of readings to be taken and the time period (in one-sixtieth of seconds) between samples. These values are part of the data file stored by SAMPLE. Subroutine UPDATE performs the command 'R'.

Command 'L'

This command lets the user set limits for an audible alarm that monitors the input voltages of the channels that are active for sampling. If the input voltage is not within the specified range, the program sounds the bell on the terminal

each second. With the correct limits, the alarm can be used to notify the user when a bridge needs rebalancing before starting an experiment. Subroutine UPDATE performs the command 'L'.

Command 'E'

This command allows the user to define the parameters for operating under automatic control. The time interval between experiments, the number of experiments to be performed and the time until the first experiment can all be set to any value. Subroutine UPDATE performs the command 'E'.

Command 'A'

This command enables the automatic operation of the program, as defined by the automatic control parameters. It is performed by the main program and subroutine AUTO. When called, subroutine AUTO will call SAMPLE and CALCULATE as needed, increment the filenames of the active channels and update the automatic operation parameters. Using the 'A' command, the operator can free himself from having to continuously monitor and initiate experiments or calculations.

Command 'B'

This command disables automatic operation of the program. The time until the next experiment will continue to count down, but AUTO will not be called if the count reaches zero.

Using the 'B' command, the user can exit automatic control, change a system parameter, reenter automatic control and not lose the correct spacing between experiments. The main program performs the 'B' command.

Command 'C'

This command performs the desired calculation on each channel that is active for calculate. Subroutine CALC performs this command but does not perform the calculations. CALC reads and prepares the data file for a subroutine called CRUNCH, which is to perform the actual calculations. CRUNCH can either be a thermal conductivity or blood perfusion routine that can be linked to CALC for use.

Command 'D'

This command allows the user to examine a data file by listing or graphing. Listing a data file consists of labeling, formatting and printing the system configuration at the time of sampling and the entire presample and sample data.

The graphics package plots the cooling curve and labels it with pertinent information from the data file. The user has the options of changing the size and position of the graph as well as expanding the time scale to display only a portion of the cooling curve.

Subroutine DATA performs the command 'D'.

Command 'S'

This command causes the computer to perform the sampling as defined by the system parameters. For each active channel the appropriate presample is collected, the probe is heated, the desired sample data is collected, and a data file is created on disk. The format of the data file is in appendix B. It consists of all possible information that might be desired at a later time to analyze the results of the experiment. Subroutine SAMP performs the command 'S'.

Command 'Z'

This command terminates the operation of the program and returns control to the RT11 monitor.

Auxiliary Program

The auxiliary program PROBEC.FOR is used to maintain the data file that contains the current probe calibration constants. The user has the option of entering calibration data for a new probe or changing calibration data for an existing probe in the file. Each time new data is entered into this file, the user can also enter the date of calibration. The probe calibration and the date are stored in this for use by the system software. Specifications for the data file PROBE.DAT are contained in appendix B.

V. RESULTS

The data acquisition system has been successfully used in a variety of configurations.

Figures 8 through 14 are actual output from the system. Figure 8 shows the cooling curve of a probe that was placed in a 100 ml solution of glycerin at 38 degrees Celsius. The presample data is shown from minus five to minus two seconds. The plot disappears from minus two to zero seconds, since that is during the heat pulse and the bridge output is not valid. The sample data is shown from zero to fifteen seconds.

This curve shows that the probe, the bridge, the signal conditioning, and the analog to digital converter all worked properly to deliver the appropriate data to the computer. It also shows that the software is capable of procuring the data, storing it, and retrieving it for display.

Figure 9 is a cooling curve from a similar experiment but the heating pulse has been reduced to only one-sixth of a second. The curve decays much more rapidly since less heat was deposited by the probe. In order to display this data in a more visible manner, the software allows for expansion of the time scale so that any portion of the curve can be shown. An example of this capability is shown in Figure 10 using the same data file as Figure 9.

Figure 11 is from a data file that was created from an experiment using a live rabbit. The blood perfusion in the




Figure 8. Thermal Conductivity Cooling Curve

GLYCERYN THERMAL CONDUCTIVITY

FILENAME: AAA002 CHANNEL: 1

PROBE: 24 HEAT PULSE: 0:10

SAMPLE TIME- 11:23: 5 DATE- 3/31/62




Figure 9. Thermal Conductivity Cooling Curve
Short Heat Pulse Duration

GLYCERYN THERMAL CONDUCTIVITY

FILENAME: AAR006 CHANNEL: 1

PROBE: 24 HEAT PULSE: 0:10

SAMPLE TIME- 11:23: 5 DATE- 3/31/82




Figure 10. Thermal Conductivity Cooling Curve

Expanded Time Scale




Figure 11. Kidney Blood Perfusion Cooling Curve

kidney could be determined from this data.

Figures 12 and 13 are from data files created from experiments conducted on the liver of a live rabbit. The probe for Figure 12 was placed very close to a large blood vessel in the liver. This results in a cooling curve that decays quite rapidly. Also present on this curve is some respiratory artifact that causes the curve to rise and fall periodically. Figure 13 is from a test on the same liver but the probe has been moved away from any blood vessels. Perfusion is still present, but at a much lower rate and the respiratory artifact is no longer visible. In order to compare these two curves, Figure 14 was generated by expanding and overlaying the data from the two files on the same plot. The difference between the high perfusion curve and the low perfusion curve is now seen as different rates of decay in the cooling curves. The higher perfusion near the large blood vessel rapidly removes heat from the vicinity of the probe, resulting in a faster decay than is seen in the low perfusion case.

RABBIT LIVER PERFUSION TEST

FILENAME: RLP001 CHANNEL: 1

PROBE: 24 HEAT PULSE: 4: 0

SAMPLE TIME- 14: 5: 7 DATE- 4/16/82




Figure 12. Rabbit Liver, High Blood Perfusion




Figure 13. Rabbit Liver, Low Blood Perfusion




Figure 14. Rabbit Liver, Comparison of High and Low Perfusion

VI. CONCLUSIONS AND RECOMMENDATIONS

Conclusions

An automated control and data acquisition system for the thermal pulse-decay method of blood perfusion measurement has been designed and tested. The system meets or exceeds the desired specifications, providing a useable tool for the continuation of research in local blood perfusion measurement.

Recommendations

Based on system performance until this point, the following recommendations are proposed as future modifications to the system. Additional desired modifications will become apparent as the system is utilized more frequently.

- 1) Incorporate into the auxiliary program, PROBEC.FOR, the capability to automatically calibrate the thermistor probes.
- 2) Enhance the data presentation, particularly the ability to show more than one cooling curve on the same plot.
- 3) Provide more disk file security to prevent accidental overwriting of data.
- 4) Provide for easy transfer of data files to or from this computer and another.
- 5) Incorporate computer routines to analyse the cooling curves for both thermal conductivity and blood perfusion.

APPENDIX A SYSTEM HARDWARE CONNECTIONS

This appendix contains the wiring list for connecting the DEC DRV11 parallel input-output board to the bridge circuit. Twotwenty-five wire ribbon cables carry the signals listed below in Table 2.

Table 2. Pin Assignments for Bridge Control

<u>Signal</u>	<u>Channel</u>	<u>Connector, Pin</u>
Presample LED	1	J1, 11
Presample LED	2	J1, 12
Presample LED	3	J1, 13
Presample LED	4	J1, 14
Presample LED	5	J1, 15
Presample LED	6	J1, 16
Heat Pulse and LED	1	J2, 11
Heat Pulse and LED	2	J2, 12
Heat Pulse and LED	3	J2, 13
Heat Pulse and LED	4	J2, 14
Heat Pulse and LED	5	J2, 15
Heat Pulse and LED	6	J2, 16
Sample LED	ALL	J1, 17

APPENDIX B
DATA FILE STRUCTURES

This appendix contains the structure of the two types of data files used by the operating software. PROBE.DAT contains the calibration data for the thermistor probes and is outlined in Table 3. The structure of the data files created by SAMPLE is outlined in Table 4.

Table 3. PROBE.DAT Structure

<u>Record</u>	<u>Variable Description</u>	<u>Disk Format</u>
First	Number of Probes	I*2
Second	Probe Number	I*2
	Month	I*2
	Day	I*2
	Year	I*2
	A	R*4
	B	R*4
	RB0	R*4
	BETA	R*4

Repeating for each probe in the data file.

Table 4. Data File Structure

<u>Variable Description</u>	<u>Variable Name</u>	<u>Disk Format</u>
Data Version	1	I*2
Time	ISTIM	I*4
Date	KDATE	3 I*2
Channel	I	I*2
Text	ITEXT	20 I*2
Active Files	IFILE(1)	3 I*2
Active Files	IFILE(2)	3 I*2
Active Files	IFILE(3)	3 I*2
Active Files	IFILE(4)	3 I*2
Active Files	IFILE(5)	3 I*2
Active Files	IFILE(6)	3 I*2
Probe Number	IPROBE	I*2
RBO	PPRBO	R*4
BETA	PPBETA	R*4
A	PPA	R*4
B	PPB	R*4
Probe Calibration Date	ICDATE	3 I*2
Bridge Balance Resistance	IPPRBB	I*2
Bridge Balance Time	IPPTBB	I*4
Bridge Balance Voltage	PPVBB	R*4
Heat Pulse Duration (Ticks)	IHEAT	I*2
Presample Readings	IPSR	I*2
Presample Period (Ticks)	IPSP	I*2
Sample Readings	ISR	I*2
Sample Period (Ticks)	ISP	I*2
Presample Data	IDATA	(IPSR) I*2
Sample Data	IDATA	(ISR) I*2

APPENDIX C

OPERATING PROGRAM

Figure 6 in Chapter IV is a flowchart of the main program. This appendix contains a listing of the FORTRAN code of the main program followed by flowcharts and code listings of the subroutines.

Figure 15 shows the flowchart of subroutine INIT (p. 46),
Figure 16 shows the flowchart of subroutine UPDATE (p.49-50),
Figure 17 shows the flowchart of subroutine SAMP (p. 57),
Figure 18 shows the flowchart of subroutine CALC (p. 63),
Figure 19 shows the flowchart of subroutine AUTO (p. 65),
Figure 20 shows the flowchart of subroutine DATA (p. 68),
Figure 21 shows the flowchart of PROBEC.FOR (p. 74), and
Figure 22 shows the flowchart of subroutine GET (77).


```

C      THIS IS THE MAIN PROGRAM FOR THERMAL PULSE
C      STUDIES.
C      IT WAS WRITTEN BY KURT L. BAUM

C      * * * COMMON VARIABLES * * *
      LOGICAL*1      HAZ(3),HAZA(5)
      INTEGER         I2(256)
      INTEGER*4       I4(12)
      REAL*4          R4(64)
      INTEGER         IDATA(1800)

      COMMON /HAZEL/HAZ,HAZA
      COMMON /VARBLE/I2,I4,R4
      COMMON /TEMPS/IDATA
C      * * * END COMMON BLOCK * * *

C      * * * BEGIN COMMON DESIGNATIONS * * *
      INTEGER IACTVS(6), IACTVC(6), IPROBE(6), IHEAT(6), IPPR88(6)
      INTEGER IFILE(6,3), ITEXT(6,20), ICDATE(6,3)
      INTEGER*4 IPPT88(6)
      REAL*4 PPR88(6), PPBETA(6), PPA(6), PPB(6), PPV88(6)
      REAL*4 CTEMP(6), RLTEMP(6)

      EQUIVALENCE      (ICOM, I2(1)),          (IACTVS(1), I2(2))
      EQUIVALENCE      (IPROBE(1), I2(8)),      (IHEAT(1), I2(14))
      EQUIVALENCE      (IPPR88(1), I2(20)),      (IPSR, I2(26))
      EQUIVALENCE      (IPSP, I2(27)),          (ISR, I2(28))
      EQUIVALENCE      (ISP, I2(29)),            (ITEXT(1,1), I2(30))
      EQUIVALENCE      (ICDATE(1,1), I2(150)),  (IACTVC(1), I2(168))
      EQUIVALENCE      (IFILE(1,1), I2(174)),    (IAUTOC, I2(192))
      EQUIVALENCE      (IEXREP, I2(193))
      EQUIVALENCE      (IPPT88(1), I4(1))
      EQUIVALENCE      (PPR88(1), R4(1)),        (PPBETA(1), R4(7))
      EQUIVALENCE      (PPA(1), R4(13)),         (PPB(1), R4(19))
      EQUIVALENCE      (PPV88(1), R4(25)),       (CTEMP(1), R4(31))
      EQUIVALENCE      (RLTEMP(1), R4(37)),      (RUALIM, R4(43))
      EQUIVALENCE      (RLALIM, R4(44)),         (AUTOTH, R4(45))
      EQUIVALENCE      (AUTOTU, R4(46)),         (EXINT, R4(47))

C      * * * END COMMON DESIGNATIONS * * *

C      * * * BEGIN LOCAL VARIABLES * * *
      LOGICAL*1 PRNTIM(9)
      INTEGER JTEMP(6)
      INTEGER ITYPE(30)
      INTEGER*4 KTIME
      INTEGER KTIM(2)
      REAL*4 CVOLT(6), TEMP

      EQUIVALENCE (KTIME, KTIM(1))
C      * * * END LOCAL VARIABLES * * *

C      * * * BEGIN CODE SEGMENT * * *
      HAZ(2)='200
      HAZ(3)='200
      HAZA(1)=27
      HAZA(2)=61
      HAZA(5)='200
      CALL INIT
      PRNTIM(9)='200

C      * * * PRINT PROGRAM IDENTIFICATION * * *
      TYPE 30
      HAZ(1)=25
      CALL PRINT(HAZ)
      HAZA(3)=40

```

```

      HAZA(4)=45
      CALL PRINT(HAZA)
      TYPE 1
1      FORMAT('THERMAL PULSE DECAY CONTROLLER')
      TYPE 30
      TYPE 2
2      FORMAT('15, WRITTEN BY KURT L. BAUM')
      ACCEPT 3
3      FORMAT(A)

C * * * PRINT SMORGASBOARD * * *
5      TYPE 30
      HAZ(1)=26
      CALL PRINT(HAZ)
      CALL IDATE(KMON,KDAY,KYEAR)
      ENCODE(3,100,PRNTIM) KMON,KDAY,KYEAR
100     FORMAT(I2,'/',I2,'/',I2)
      CALL PRINT('TIME')
      HAZA(3)=32
      HAZA(4)=52
      CALL PRINT(HAZA)
      CALL PRINT('DATE')
      HAZA(4)=57
      CALL PRINT(HAZA)
      CALL PRINT(PRNTIM)
      HAZA(4)=77
      CALL PRINT(HAZA)
      CALL PRINT('DATA VERSION ONE')
      TYPE 105,1,2,3,4,5,6
105     FORMAT(' CHANNEL',1,6(7X,I1))
      DO 107 I=1,6
      JTEMP(I)=IACTVS(I)
      IF(IACTVS(I).NE.0) GOTO 107
      JTEMP(I)=0
107     CONTINUE
      TYPE 110,JTEMP
110     FORMAT('O SAMPLE STATUS',1,6(6X,A2))
      DO 112 I=1,6
      JTEMP(I)=IACTVC(I)
      IF(IACTVC(I).NE.0) GOTO 112
      JTEMP(I)=0
112     CONTINUE
      TYPE 115,JTEMP
115     FORMAT(' CALCULATE STATUS',1,6(6X,A2))
      TYPE 120,((IFILE(I,K),K=1,3),I=1,6)
120     FORMAT('F FILENAME',1,6(2X,3A2))
      TYPE 135,Iprobe
135     FORMAT('P PROBE',1,6I8)
      TYPE 136
136     FORMAT(' CURRENT MILLIVOLTS')
      TYPE 137
137     FORMAT(' CURRENT TEMPERATURE')
      TYPE 138,RLTEMP
138     FORMAT(' LAST TEMP',1,6F8.3)
      TYPE 140,(IHEAT(I)/60,IHEAT(I)-IHEAT(I)/60*60,I=1,6)
140     FORMAT('H HEAT (SEC:TICKS)',1,6(15,' ',I2))
      TYPE 142
142     FORMAT(' DURATION (SEC:TIC) FREQUENCY (HZ)',
1          ' READINGS PERIOD')
      KTIM(1)=0
      KTIM(2)=IPSR*IPSP
      CALL CVTTIM(KTIME,KHRS,KMIN,KSEC,KTICKS)
      TYPE 145,KSEC,KTICKS,60/IPSP,IPSR,IPSP
145     FORMAT('R PRESAMPLE',1,12,' ',I2,' ',I4,' ',I2)
1          ' ',I2,' ',I4,' ',I2)

```

```

      KTIM(2)=ISR*ISP
      CALL CVTTIM(KTIME,KHRS,KMIN,KSEC,KTICKS)
      TYPE 150,KSEC,KTICKS,60/ISP,ISR,ISP
150  FORMAT('  SAMPLE      ',I2,' ',I2,'
      1      I2,'      ',I4,'      ',I2)
      TYPE 30
      TYPE 151,RUALIM
151  FORMAT('L  UPPER ALARM LIMIT (MV)      ',F5.2,
      1      T45,'A  AUTOMATIC RUN')
      TYPE 152,RLALIM
152  FORMAT('  LOWER ALARM LIMIT (MV)      ',F5.2,
      1      T45,'B  BREAK FROM AUTO')
      KTIM(1)=INT(EXINT/65536.0)
      KTIM(2)=INT(EXINT-65536.0*KTIM(1))
      CALL CVTTIM(KTIME,KHRS,KMIN,KSEC,KTICKS)
      TYPE 153,KMIN,KSEC
153  FORMAT('E  EXPERIMENT INTERVAL (MIN:SEC) ',I2,' ',I2,
      1      T45,'C  CALCULATE')
      TYPE 160,IEXP
160  FORMAT('  EXPERIMENT REPETITIONS      ',I4,
      1      T45,'D  LOOK AT DATA FILE')
      TYPE 163
163  FORMAT('  NEXT EXPERIMENT (MIN:SEC)      0: 0',
      1      T45,'S  SAMPLE')
      TYPE 166
166  FORMAT(T45,'Z  QUIT')
      HAZA(3)=52
      HAZA(4)=32
      CALL PRINT(HAZA)
      TYPE 190
190  FORMAT(' PLEASE ENTER YOUR COMMAND',S)
30  FORMAT(' ')
      IF(ICOM.NE.63) GOTO 50
      TYPE 30
      CALL PRINT(HAZA)
      TYPE 40
40  FORMAT('UNDER AUTOMATIC OPERATION')

C  * * * CHECK USER INPUT * * *
50  CALL IPOKE('44',"010100.OR.IPEEK('44'))
      ITEMP=ITINR()
      CALL IPOKE('44',"167677.AND.IPEEK('44'))
      IF (ITEMP.GE.0) GOTO 200

C  * * * CHECK FOR JUMP TO AUTO * * *
      IF(ICOM.EQ.63.AND.AUTOTN.EQ.0.0) GOTO 207

C  * * * CHECK FOR TIME AND TEMP UPDATE * * *
77  CALL GTIM(KTIME)
      CALL CVTTIM(KTIME,KHRS,KMIN,KSEC,KTICKS)
      IF (KSEC.EQ.ISEC) GOTO 50
      ISEC=KSEC

C  * * * REPRINT TIMES AND TEMPS * * *
80  ENCODE(8,80,PRNTIM) KHRS,KMIN,KSEC
      FORMAT(I2,' ',I2,' ',I2)

C  * * * GET CURRENT TEMPS * * *
      K=1
      DO 90 I=1,6
      CALL IPOKE('176770,K)
      J=IPEEK('176772)
      IF(J.GT.'3777) J=J-'10000
      CVOLT(I)=J/389.905
C  * * * 389 905='3777/5.25 * * *

```

```

      IF(IACTVS(I).EQ.0) GOTO 88
      IF(CVOLT(I) GE. RLALIM AND CVOLT(I) LE. RUALIM) GOTO 88
      HAZ(1)=7
      CALL PRINT(HAZ)
88      RBB=FLOAT(IPRBB(I))/2.0+1000.0
      VRS=119756.88
      RS=22075.0
      TEMP=VRS/(VRS/(RBB+RS)+CVOLT(I)/1000)-RS
C* * * WHERE 119756.88=VB*22075, VB=5.425 VOLTS * * *
      CTEMP(I)=PPA(I)-PPB(I)*ALOG(RB)
      K=K+400
90      CONTINUE

C * * * PRINT TIME * * *
      TYPE 30
      HAZA(3)=32
      HAZA(4)=38
      CALL PRINT(HAZA)
      CALL PRINT(PRNTIM)

C * * * PRINT CURRENT MILLIVOLTS * * *
      TYPE 30
      HAZA(3)=38
      HAZA(4)=56
      CALL PRINT(HAZA)
      TYPE 95,CVOLT
95      FORMAT(6F8.3)

C * * * PRINT CURRENT TEMPS * * *
      TYPE 30
      HAZA(3)=39
      CALL PRINT(HAZA)
      TYPE 97,CTEMP
97      FORMAT(6F8.3)

C * * * UPDATE AND PRINT NEXT AUTO TIME * * *
      IF(AUTOTN.EQ.0.0) GOTO 99
      IF(IAUTOC.EQ.0) GOTO 96
      TIME=65536.0*KTIM(1)+KTIM(2)
      IF(AUTOTU.GT.TIME) AUTOTU=AUTOTU-5384000.0
      AUTOTN=AUTOTN-(TIME-AUTOTU)
      IF(AUTOTN.LT.0.0) AUTOTN=0.0
      AUTOTU=TIME
96      TYPE 30
      HAZA(3)=50
      HAZA(4)=66
      CALL PRINT(HAZA)
      KTIM(1)=INT(AUTOTN/65536.0)
      KTIM(2)=INT(AUTOTN-65536.0*KTIM(1))
      CALL CVTTIM(KTIME,KHRS,KMIN,KSEC,KTICKS)
      TYPE 98,KMIN,KSEC
98      FORMAT(I2,' ',I2)

C * * * RESTORE CURSOR * * *
99      TYPE 30
      HAZA(3)=53
      HAZA(4)=59
      CALL PRINT(HAZA)
      GOTO 50

C * * * PERFORM RECIEVED COMMAND * * *
200      ICOM=ITEMP
205      TYPE 30
207      IF (ICOM.NE.65) GOTO 215
      IF (IAUTOC.NE.0) GOTO 210

```

```
IAUTOC=1
CALL GTIM(KTIME)
AUTOTU=65536 0:KTIME(1)+KTIME(2)
210 IF(AUTOTU.EQ.0) CALL AUTO
215 IF (ICOM.EQ.67) CALL CALC
IF (ICOM.EQ.68) CALL DATA
IF (ICOM.EQ.90) GOTO 300
IF (ICOM.EQ.93) CALL SAMP
220 IF (ICOM.GT.68.AND.ICOM.LT.83) CALL UPDATE
GOTO 5
300 STOP '      END    OF    PROGRAM    '
END
```




Figure 15. Flowchart of Subroutine INIT

```

SUBROUTINE INIT

C      THIS SUBROUTINE INITIALIZES THE SYSTEM PARAMETERS

C * * * COMMON VARIABLES * * *
LOGICAL*1      HAZ(3),HAZA(5)
INTEGER        I2(256)
INTEGER*4      I4(12)
REAL*4         R4(64)
INTEGER        IDATA(1800)

COMMON /HAZEL/HAZ,HAZA
COMMON /VARBLE/I2,I4,R4
COMMON /TEMPS/IDATA
C * * * END COMMON BLOCK * * *

C * * * BEGIN COMMON DESIGNATIONS * * *
INTEGER IACTVS(6),IACTVC(6),IPROBE(6),IHEAT(6),IPRBB(6)
INTEGER IFILE(6,3),ITEXT(6,20),ICDATE(6,3)

EQUIVALENCE      (ICOM,I2(1)),          (IACTVS(1),I2(2))
EQUIVALENCE      (IPROBE(1),I2(8)),      (IHEAT(1),I2(14))
EQUIVALENCE      (IPRBB(1),I2(20)),      (IPSR,I2(26))
EQUIVALENCE      (IPSP,I2(27)),          (ISR,I2(29))
EQUIVALENCE      (ISP,I2(29)),            (ITEXT(1,1),I2(30))
EQUIVALENCE      (ICDATE(1,1),I2(150)),  (IACTVC(1),I2(163))
EQUIVALENCE      (IFILE(1,1),I2(174)),   (IAUTOC,I2(193))
EQUIVALENCE      (IEXREP,I2(193))
EQUIVALENCE      (PPYBB(1),R4(25)),       (CTEMP(1),R4(31))
EQUIVALENCE      (RLTEMP(1),R4(37)),      (RUALIM,R4(43))
EQUIVALENCE      (RLALIM,R4(44)),         (AUTOTH,R4(45))
EQUIVALENCE      (AUTOTU,R4(46)),         (EXINT,R4(47))

C * * * END COMON DESIGNATIONS * * *

C * * * BEGIN CODE SEGMENT * * *
DO 100 I=1,256
  I2(I)=0
100 CONTINUE

DO 110 I=1,64
  R4(I)=0.0
110 CONTINUE

DO 120 I=1,12
  I4(I)=0
120 CONTINUE

DO 5 I=1,6
  IFILE(1,1)='AA'
  IFILE(1,3)='00'
  CONTINUE
  IFILE(1,2)='A0'
  IFILE(2,2)='B0'
  IFILE(3,2)='C0'
  IFILE(4,2)='D0'
  IFILE(5,2)='E0'
  IFILE(6,2)='F0'
5 CONTINUE

```

```
DO 10 I=1.6  
  IACTVSC(I)= 'ON'  
  IACT'CCI)= 'ON'  
  IHEAT(I)=120  
10  CONTINUE  
  
DO 20 J=1.6  
DO 20 I=1.20  
  ITEXT(J,I)= '52124'  
20  CONTINUE  
  
  IPSR=10  
  IPSP=60  
  ISR=10  
  ISP=60  
  
  RUALIM=6.0  
  RLALIM=-6.0  
  
  EXINT=60.0  
  IEXREP=1  
  AUTOTN=30.0  
  
  RETURN  
  END
```





Figure 16. Flowchart of Subroutine UPDATE




Figure 16. Continued

```

SUBROUTINE UPDATE

C      THIS SUBROUTINE ALLOWS THE USER TO CHANGE
C      PARAMETERS TO TAILOR SYSTEM CONFIGURATION

C      * * * COMMON VARIABLES * * *
      LOGICAL*1      HAZ(3),HAZA(5)
      INTEGER        I2(256)
      INTEGER*4      I4(12)
      REAL*4         R4(64)
      INTEGER        IDATA(1800)

      COMMON /HAZEL/HAZ,HAZA
      COMMON /VARBLE/I2,I4,R4
      COMMON /TEMPS/IDATA
C      * * * END COMMON BLOCK * * *

C      * * * BEGIN COMMON DESIGNATIONS * * *
      INTEGER IACTVS(6), IACTVC(6), IPROBE(6), IHEAT(6), IPPR88(6)
      INTEGER IFILE(6,3), ITEXT(6,20), ICDATE(6,3)
      INTEGER*4 IPPT88(6)
      REAL*4 PPR88(6), PPBETA(6), PPA(6), PPB(6), PPVBB(6)
      REAL*4 CTEMP(6), RLTEMP(6)

      EQUIVALENCE      (ICOM, I2(1)),          (IACTVS(1), I2(2))
      EQUIVALENCE      (IPROBE(1), I2(8)),      (IHEAT(1), I2(14))
      EQUIVALENCE      (IPPR88(1), I2(28)),      (IPSR, I2(26))
      EQUIVALENCE      (IPSP, I2(27)),          (ISR, I2(28))
      EQUIVALENCE      (ISP, I2(29)),           (ITEXT(1,1), I2(30))
      EQUIVALENCE      (ICDATE(1,1), I2(150)),   (IACTVC(1), I2(168))
      EQUIVALENCE      (IFILE(1,1), I2(174)),    (IAUTOC, I2(192))
      EQUIVALENCE      (IEXREP, I2(193))
      EQUIVALENCE      (IPPT88(1), I4(1))
      EQUIVALENCE      (PPR88(1), R4(1)),        (PPBETA(1), R4(7))
      EQUIVALENCE      (PPA(1), R4(13)),        (PPB(1), R4(19))
      EQUIVALENCE      (PPVBB(1), R4(25)),      (CTEMP(1), R4(31))
      EQUIVALENCE      (RLTEMP(1), R4(37)),     (RUALIM, R4(43))
      EQUIVALENCE      (RLALIM, R4(44)),        (AUTOTN, R4(45))
      EQUIVALENCE      (AUTOTU, R4(46)),        (EXINT, R4(47))
C      * * * END COMMON DESIGNATIONS * * *

C      * * * BEGIN LOCAL VARIABLES * * *
      INTEGER JDATA(12)
      EQUIVALENCE(ATEM, JDATA(5)), (STEM, JDATA(7))
      EQUIVALENCE(RB0TEM, JDATA(9)), (BETAT, JDATA(11))
      EQUIVALENCE(IYEAR, JDATA(4))
      EQUIVALENCE(IMONTH, JDATA(2)), (IDAY, JDATA(3))
C      * * * END LOCAL VARIABLES * * *

C      * * * BEGIN CODE SEGMENT * * *

C      * * * CLEAR ROW AND GET SET TO TYPE REQUEST * * *
      HAZA(3)=52
      HAZA(4)=32
      CALL PRINT(HAZA)
      HAZ(1)=27
      HAZ(2)=89
      CALL PRINT(HAZ)
20  FORMAT( )

C      * * * CHANGE HEAT PULSE DURATIONS * * *
      IF (ICOM.NE 72) GOTO 100

```

```

      TYPE 30
30      FORMAT('UPDATE HEAT PULSE DURATION OF WHICH CHANNEL ? ', $)
      ACCEPT 43, I
40      FORMAT( I )
      TYPE 20
      IF ( I.LT.1.OR.I.GT.6 ) GOTO 5
      CALL PRINT(HAZA)
      CALL PRINT(HAZ)
      TYPE 50
33      FORMAT('ENTER NEW HEAT PULSE DURATION ', $)
      ACCEPT 60, IHEAT(I)
60      FORMAT( I3 )
      GOTO 1000

```

```

C * * * CHANGE FILENAME * * *
100      IF (ICOM.NE.70) GOTO 200
      TYPE 110
110      FORMAT('UPDATE FILENAME OF WHICH CHANNEL ? ', $)
      ACCEPT 120, I
120      FORMAT( I )
      TYPE 20
      IF ( I.LT.1.OR.I.GT.6 ) GOTO 5
      CALL PRINT(HAZA)
      CALL PRINT(HAZ)
      TYPE 130
130      FORMAT('ENTER NEW FILENAME - - - AAAAAA ', $)
140      ACCEPT 150, IFILE(I,1), IFILE(I,2), IFILE(I,3)
150      FORMAT(3A2)
      TYPE 20
      DECODE(2, 160, IFILE(I,3), ERR=170) IFILEA
160      FORMAT( I2 )
      IF ( IFILE(I,2).GT."30000".AND.IFILE(I,2).LT."35000" ) GOTO 1000
170      CALL PRINT(HAZA)
      CALL PRINT(HAZ)
      TYPE 190
190      FORMAT('INCORRECT ENTRY - TRY AGAIN - AAAAAA ', $)
      GOTO 140

```

```

C * * * CHANGE PROBES, TEXTS, AND BRIDGE BALANCES * * *
200      IF (ICOM.NE.90) GOTO 400
      HAZ(1)=26
      HAZ(2)="200"
205      CALL PRINT(HAZ)
      TYPE 20
      TYPE 210
210      FORMAT('UPDATE PROBE INFORMATION OF WHICH CHANNEL ? ', $)
      ACCEPT 215, I
215      FORMAT( I )
      TYPE 20
      IF ( I.LT.1.OR.I.GT.6 ) GOTO 205
217      CALL PRINT(HAZ)
      TYPE 210, I
218      FORMAT('PROBE IS ON CHANNEL - ', I)
      TYPE 220, IPROBE(I)
220      FORMAT('PROBE NUMBER          - ', I3)
      TYPE 225, PPRB0(I)
225      FORMAT('RB0                      - ', F8.2)
      TYPE 230, PPBETA(I)
230      FORMAT('BETA                      - ', F8.3)
      TYPE 235, PPA(I)
235      FORMAT('A                        - ', F8.4)

```

```

240 TYPE 240,PPB(I)
   FORMAT('B' - ',FS,4)
   TYPE 245,(ICDATE(I,K),K=1,3)
245 FORMAT('PROBE CALIBRATED ON - ',I2,',',I2,',',I2)
   TYPE 250,IPRBB(I)
250 FORMAT('BRIDGE BALANCE RESISTANCE - - ',I4)
   TYPE 255,PPVB(I)
255 FORMAT('BRIDGE BALANCE VOLTAGE - - ',F5.2)
   CALL CVTTIM(IPPTBB(I),KHRS,KMIN,KSEC,KTICKS)
   TYPE 260,KHRS,KMIN,KSEC
260 FORMAT('BRIDGE BALANCED AT - - ',I2,',',I2,',',I2)
   TYPE 265,(ITEXT(I,K),K=1,20)
265 FORMAT('TEXT - - ',20A2)
   TYPE 20
   TYPE 20
   TYPE 270
270 FORMAT('TYPE "P" TO CHANGE PROBE')
   TYPE 275
275 FORMAT('      "B" TO CHANGE BALANCE')
   TYPE 280
280 FORMAT('      "T" TO CHANGE TEXT')
   TYPE 283
283 FORMAT('      "C" TO WORK WITH NEW CHANNEL')
   TYPE 285
285 FORMAT('      OR "E" TO RETURN TO SMORGASBOARD ',)
288 ACCEPT 290,K
290 FORMAT('A ')
   IF (K.EQ.'E') GOTO 1000
C * * * CHANGE PROBE NUMBER AND CALIBRATION DATA * * *
   IF (K.NE.'P') GOTO 330
   TYPE 295
295 FORMAT('ENTER NEW PROBE NUMBER ',)
   CALL ASSIGN(3,'F00:PROBE.DAT',13)
   READ (3)JNUM
   INUM=JNUM*12
   READ(3)(IDATA(K),K=1,INUM)
   CALL CLOSE(3)
   ACCEPT 290,K
290 FORMAT('I3)
   J=1
300 IF(K.EQ.IDATA(J)) GOTO 310
   J=J+12
   IF(J.LT.INUM) GOTO 300
   TYPE 305
305 FORMAT('NO PROBE EXISTS WITH THAT DESIGNATION')
   GOTO 217
310 IPRBB(I)=K
   DO 320 K=1,12
   JOATA(K)=IDATA(J)
   J=J+1
320 CONTINUE
   PPAC(I)=ATEM
   PPB(I)=BTEM
   PPRBB(I)=RBBTEM
   PPSETA(I)=BETAT
   ICDATE(I,1)=IMONTH
   ICDATE(I,2)=IDAY
   ICDATE(I,3)=IYEAR
   GOTO 217
C * * * CHANGE BRIDGE BALANCE * * *
330 IF (K.NE.'B') GOTO 370
   TYPE 335
335 FORMAT('CURRENT RAW DATA          CURRENT MILLIVOLTS')
   K=1
   IF(I.EQ.1) GOTO 345

```

```

DO 340 J=1,I-1
K=K+*400
340 CONTINUE
345 HAZA(3)=51
HAZA(4)=32
347 TYPE 20
CALL PRINT(HAZA)
ISEC=KSEC
CALL IPOKE('176770,K)
J=IPEEK('176772)
J1=J
IF(J.GT.'3777) J1=J-'10000
PPV88(I)=J1/389.9047619
TYPE 350,J,PPV88(I)
350 FORMAT(SX.06,10X,F7.4)
TYPE 355
355 FORMAT('PRESS ANY KEY WHEN DESIRED BALANCE IS REACHED ')
360 CALL IPOKE('44,'010100.OR.IPEEK('44))
J=ITINR()
CALL IPOKE('44,'167677.AND.IPEEK('44))
IF(J.GE.0) GOTO 365
CALL GTIM(IPPT88(I))
CALL CVTTIM(IPPT88(I),KHRS,KMIN,KSEC,KTICKS)
IF(KSEC.EQ.ISEC) GOTO 360
GOTO 347
365 TYPE 366
366 FORMAT('ENTER BRIDGE RESISTANCE BALANCE ',*)
ACCEPT 367,IPPR88(I)
367 FORMAT(I4)
GOTO 217
C * * * CHANGE TEXT * * *
370 IF(K.NE.'T') GOTO 385
TYPE 375
375 FORMAT('ENTER NEW TEXT - 12345678901234567890',
1 '12345678901234567890')
TYPE 380
380 FORMAT(' ',*)
ACCEPT 382,(ITEXT(I,K),K=1,20)
382 FORMAT(20A2)
GOTO 217
385 IF(K.EQ.'C') GOTO 285
TYPE 390
390 FORMAT('INVALID ENTRY - - TRY AGAIN',*)
GOTO 280

C * * * CHANGE SAMPLE PARAMETERS * * *
400 IF(ICOM.NE.92) GOTO 500
TYPE 410
410 FORMAT('ENTER NEW VALUE OF PRESAMPLE DURATION (SEC) ',*)
ACCEPT 420,IPSR
420 FORMAT(I4)
TYPE 20
CALL PRINT(HAZA)
CALL PRINT(HAZ)
TYPE 430
430 FORMAT('ENTER NEW PRESAMPLE FREQUENCY (HZ) ',*)
ACCEPT 420,IPSP
IPSP=60/IPSP
IPSR=IPSR*60/IPSP
TYPE 20
CALL PRINT(HAZA)
CALL PRINT(HAZ)
TYPE 440

```

```

140  FORMAT('ENTER NEW VALUE OF SAMPLE DURATION (SEC) ',%)
      ACCEPT 420,ISR
      TYPE 20
      CALL PRINT(HAZA)
      CALL PRINT(HAZ)
      TYPE 450
150  FORMAT('ENTER NEW SAMPLE FREQUENCY (HZ) ',%)
      ACCEPT 420,ISP
      ISP=60./ISP
      ISR=ISR*59/ISP
      GOTO 1000

C * * * CHANGE CHANNEL ACTIVITIES * * *
500  IF(ICOM.NE.79) GOTO 600
      TYPE 510
510  FORMAT('ENTER SAMPLE CHANNEL ACTIVITY AS A SIX BIT BINARY',
1      ' NUMBER ',%)
      ACCEPT 520,IACTVS
520  FORMAT(601)
      TYPE 20
      CALL PRINT(HAZA)
      CALL PRINT(HAZ)
      TYPE 540
540  FORMAT('ENTER CALCULATE CHANNEL ACTIVITY AS A SIX BIT',
1      ' BINARY NUMBER ',%)
      ACCEPT 550,IACTVC
550  FORMAT(601)
      CC 570 J=1,6
      IF(IACTVS(J).EQ.0) GOTO 560
      IACTVS(J)='ON'
560  IF(IACTVC(J).EQ.0) GOTO 570
      IACTVC(J)='ON'
570  CONTINUE

C * * * CHANGE ALARM LIMITS * * *
500  IF(ICOM.NE.76) GOTO 700
      TYPE 610
610  FORMAT('ENTER NEW UPPER ALARM LIMIT (MV) ',%)
      ACCEPT 620,RUALIM
520  FORMAT(F5.2)
      TYPE 20
      CALL PRINT(HAZA)
      CALL PRINT(HAZ)
      TYPE 630
630  FORMAT('ENTER NEW LOWER ALARM LIMIT (MV) ',%)
      ACCEPT 640,RLALIM
540  FORMAT(F5.2)
      GOTO 1000

C * * * CHANGE EXPERIMENT REPEAT DATA * * *
700  IF(ICOM.NE.69) GOTO 800
      TYPE 710
710  FORMAT('ENTER NEW EXPERIMENT INTERVAL (SEC) ',%)
      ACCEPT 720,J
720  FORMAT(I4)
      EXINT=J*60.0
      TYPE 20
      CALL PRINT(HAZA)
      CALL PRINT(HAZ)

```

```

      TYPE 730
      FORMAT('ENTER NEW EXPERIMENT REPETITIONS ',1)
      ACCEPT 740,IEKREP
      740  FORMAT(I3)
           TYPE 20
           CALL PRINT(HAZA)
           CALL PRINT(HAZ)
           TYPE 750
      750  FORMAT('ENTER NEW TIME UNTL NEXT EXPERIMENT (SEC) ',1)
           ACCEPT 760,J
      760  FORMAT(I4)
           AUTOTN=J*60
           IAUTO=0
           GOTO 1000

      300  CONTINUE
      1000 HAZ(2)="200
           RETURN
           END
```





Figure 17. Flowchart of Subroutine SAMP

```

SUBROUTINE SAMP

C      THIS SUBROUTINE PERFORMS THE ACTUAL SAMPLING OF THE
C      DATA AND CREATES DATA FILES ON DISK FOR EACH ACTIVE
C      CHANNEL

C      * * * COMMON VARIABLES * * *
      LOGICAL*1      HAZ(3),HAZA(5)
      INTEGER         I2(256)
      INTEGER*4       I4(12)
      REAL*4          R4(64)
      INTEGER         IDATA(1800)

      COMMON /HAZEL/HAZ,HAZA
      COMMON /VARBLE/I2,I4,R4
      COMMON /TEMPS/IDATA
C      * * * END COMMON BLOCK * * *

C      * * * BEGIN COMMON DESIGNATIONS * * *
      INTEGER IACTVS(6),IACTVC(6),IPROBE(6),IHEAT(5),IPPRB8(6)
      INTEGER IFILE(6,3),ITEXT(6,20),ICDATE(6,3)
      INTEGER*4 IPPTB8(6)
      REAL*4 PPRB8(6),PPBETA(6),PPA(6),PPB(6),PPV8(6)
      REAL*4 CTEMP(6),RLTEMP(6)

      EQUIVALENCE      (ICOM,I2(1)),          (IACTVS(1),I2(2))
      EQUIVALENCE      (IPROBE(1),I2(8)),      (IHEAT(1),I2(14))
      EQUIVALENCE      (IPPRB8(1),I2(20)),      (IPSR,I2(26))
      EQUIVALENCE      (IPSP,I2(27)),          (ISR,I2(28))
      EQUIVALENCE      (ISP,I2(29)),            (ITEXT(1,1),I2(30))
      EQUIVALENCE      (ICDATE(1,1),I2(150)),   (IACTVC(1),I2(168))
      EQUIVALENCE      (IFILE(1,1),I2(174)),     (IAUTOC,I2(192))
      EQUIVALENCE      (IPPTB8(1),I4(1))
      EQUIVALENCE      (PPRB8(1),R4(1)),        (PPBETA(1),R4(7))
      EQUIVALENCE      (PPA(1),R4(13)),         (PPB(1),R4(19))
      EQUIVALENCE      (PPV8(1),R4(25)),        (CTEMP(1),R4(31))
      EQUIVALENCE      (RLTEMP(1),R4(37)),      (RUALIM,R4(43))

C      * * * END COMON DESIGNATIONS * * *

C      * * * BEGIN LOCAL DECLARATIONS * * *
      INTEGER IOFFST(6)
      INTEGER ITIM(6)
      INTEGER*2 KTIM(2)
      INTEGER ITIM8(6)
      INTEGER IFILET(18)
      INTEGER IFILEA(7)
      INTEGER KDATE(3)
      INTEGER*4 JTIM,ISTIM
      EQUIVALENCE (KMON,KDATE(1)),(KDAY,KDATE(2))
      EQUIVALENCE (KYEAR,KDATE(3)),(JTIM,KTIM(1))
C      * * * END LOCAL DECLARATIONS * * *

C      * * * BEGIN CODE SEGMENT * * *
      HAZ(1)=26
      TYPE 1
      FORMAT('NOW PERFORMING SAMPLE')
      CALL PRINT(HAZ)

C      * * * CHANGE LAST TEMPERATURES * * *
      DO 10 I=1,6
      RLTEMP(I)=CTEMP(I)
10      CONTINUE

C      * * * CLEAR IDATA * * *
      DO 5 I=1,1800

```

```

      IDATA(I)=1
5      CONTINUE

C *** SET POINTERS INTO IDATA ***
      JTEMP=2+IPSR+ISR
      K=1
      DO 20 I=1,6
      IOFFST(I)=IPSR+2
      IF (IACTVS(I).EQ.0) GOTO 20
      IOFFST(I)=K
      K=K+JTEMP
20      CONTINUE
      J1=IOFFST(1)
      J2=IOFFST(2)
      J3=IOFFST(3)
      J4=IOFFST(4)
      J5=IOFFST(5)
      J6=IOFFST(6)

C *** FIGURE OUT SHORTEST AND LONGEST HEATING PULSES ***
      IHEATS=1200
      DO 60 I=1,6
      IF (IACTVS(I).EQ.0) GOTO 60
      IF (IHEATS.LE.IHEAT(I)) GOTO 60
      IHEATS=IHEAT(I)
60      CONTINUE
      IHEATL=0
      DO 70 I=1,6
      IF (IACTVS(I).EQ.0) GOTO 70
      IF (IHEATL.GE.IHEAT(I)) GOTO 70
      IHEATL=IHEAT(I)
70      CONTINUE

C *** SET TIMES TO PERFORM PRESAMPLE ***
      ITIMA=1
      ITIME=IPSR*IPSP+ITIMA+IHEATL
      ITIMO=ITIME-IHEATS
      DO 75 I=1,6
      ITIMC(I)=ITIME
      ITIMB(I)=ITIME
75      CONTINUE
      DO 80 I=1,6
      IF (IACTVS(I).EQ.0) GOTO 80
      ITIMC(I)=ITIME-IHEAT(I)
      ITIMB(I)=ITIMC(I)-IPSR*IPSP
80      CONTINUE
      ITIME=ITIME+ISR*ISP
      ICOUNT=1

C *** CLEAR OUTPUT VARIABLE ***
      IPID=0

C *** SET TIME OFFSET ***
      CALL GTIM(JTIM)
      ITOFF=KTIM(2)

C *** PERFORM EXPERIMENT ***
90      CALL GTIM(JTIM)
      ITIM=KTIM(2)-ITOFF
      IF (ITIM.NE.ITIMA) GOTO 90
      ITIMA=ITIMA+1
      CALL IPOKE("176770","1")
      IDATA(J1)=IPEEK("176772")
      CALL IPOKE("176770","101")
      IDATA(J2)=IPEEK("176772")

```

```

CALL IPOKE("176770","1001)
IDATA(J3)=IPEEK("176772)
CALL IPOKE("176770","1401)
IDATA(J4)=IPEEK("176772)
CALL IPOKE("176770","2001)
IDATA(J5)=IPEEK("176772)
CALL IPOKE("176770","2401)
IDATA(J6)=IPEEK("176772)

IF(ITIM.GE.ITIME) GOTO 220

IF(ITIM.NE.ITIMB(1)) GOTO 110
IF(ITIM.EQ.ITIMC(1)) GOTO 100
IPIO=IPIO.OR."1
J1=J1+1
ITIMB(1)=ITIMB(1)+IPSP
GOTO 110
100 IPIO="177776.AND.IPIO.OR."400
110 IF(ITIM.NE.ITIMB(2)) GOTO 130
IF(ITIM.EQ.ITIMC(2)) GOTO 120
IPIO=IPIO.OR."2
J2=J2+1
ITIMB(2)=ITIMB(2)+IPSP
GOTO 130
120 IPIO="177775.AND.IPIO.OR."1000
130 IF(ITIM.NE.ITIMB(3)) GOTO 150
IF(ITIM.EQ.ITIMC(3)) GOTO 140
IPIO=IPIO.OR."4
J3=J3+1
ITIMB(3)=ITIMB(3)+IPSP
GOTO 150
140 IPIO="177773.AND.IPIO.OR."2000
150 IF(ITIM.NE.ITIMB(4)) GOTO 170
IF(ITIM.EQ.ITIMC(4)) GOTO 160
IPIO=IPIO.OR."10
J4=J4+1
ITIMB(4)=ITIMB(4)+IPSP
GOTO 170
160 IPIO="177767.AND.IPIO.OR."4000
170 IF(ITIM.NE.ITIMB(5)) GOTO 190
IF(ITIM.EQ.ITIMC(5)) GOTO 180
IPIO=IPIO.OR."20
J5=J5+1
ITIMB(5)=ITIMB(5)+IPSP
GOTO 190
180 IPIO="177757.AND.IPIO.OR."10000
190 IF(ITIM.NE.ITIMB(6)) GOTO 210
IF(ITIM.EQ.ITIMC(6)) GOTO 200
IPIO=IPIO.OR."40
J6=J6+1
ITIMB(6)=ITIMB(6)+IPSP
GOTO 210
200 IPIO="177737.AND.IPIO.OR."20000
210 CALL IPOKE("177556,IPIO)
GOTO 90
220 CALL IPOKE("177556,"100)
IF(ITIM.EQ.ITIME) GOTO 300
IF(ITIM.EQ.ITIME) CALL GTIM(ISTIM)

```

```

        ICOUNT=ICOUNT-1
        IF(ICOUNT.NE 0) GOTO 90
        ICOUNT=ISP
        J1=J1+1
        J2=J2+1
        J3=J3+1
        J4=J4+1
        J5=J5+1
        J6=J6+1
        GOTO 90

C * * * CLEAR OUTPUT LIGHTS * * *
300      CALL IPOKE('177556.0')

C * * * OUTPUT TEMPORARY PRESAMPLE DATA * * *
        TYPE 310
310      FORMAT('EXPERIMENT IS DONE')
        TYPE 320
320      FORMAT('FIRST TEN PRESAMPLE POINTS OF EACH CHANNEL')
        DO 345 I=1,6
        TYPE 340, (IDATA(J), J=IOFFST(I), IOFFST(I)+9)
340      FORMAT(1007)
345      CONTINUE
        TYPE 410
410      FORMAT( )
        TYPE 420, (IOFFST(I), I=1,6)
420      FORMAT(618)

C * * * OUTPUT TEMPORARY SAMPLE DATA * * *
        TYPE 410
        TYPE 550
550      FORMAT('FIRST TEN SAMPLE POINTS OF EACH CHANNEL')
        TYPE 410
        DO 570 I=1,6
        TYPE 560, (IDATA(J), J=IOFFST(I)+IPSR, IOFFST(I)+IPSR+9)
560      FORMAT(1007)
570      CONTINUE

C * * * OUTPUT DATA FILES * * *
        IFILEA(1)='FD'
        IFILEA(2)='1.'
        IFILEA(6)='D'
        IFILEA(7)='AT'
        CALL IDATE(KMON, KDAY, KYEAR)
        DO 590 I=1,18
        IFILET(I)='20040'
590      CONTINUE
        DO 610 I=1,6
        IF (IACTVS(I).EQ.0) GOTO 610
        DO 600 J=1,3
        K=I*3-3+J
        IFILET(K)=IFILE(I,J)
600      CONTINUE
610      CONTINUE
        DO 700 I=1,6
        IF (IACTVS(I).EQ.0) GOTO 700
        DO 620 J=1,3
        IFILEA(J+2)=IFILE(I,J)
620      CONTINUE
        CALL ASSIGN(3, IFILEA, 14)
        WRITE(3) 1, ISTIM, KDATE, I, IFILET, (ITEXT(I,K), K=1,20),
1          IPROBE(I), PPR88(I), PPBETA(I), PPA(I), PPS(I),
2          (ICDATE(I,K), K=1,3),
3          IPPR88(I), IPPT88(I), PPVBB(I), IHEAT(I), IPSR, IPSP,
          ISR, ISP, (IDATA(K), K=IOFFST(I), IOFFST(I)+IPSR-ISR-1)

```

```
      CALL CLOSE(3)
700    CONTINUE

C   * * *   WAIT FOR RESPONSE BEFORE RETURNING TO MAIN   * * *
      IF( ICOM.EQ.65 ) GOTO 820
      TYPE 800
800    FORMAT( 'PRESS RETURN TO EXIT BACK TO SMORGASBORD' )
      ACCEPT 810,I
810    FORMAT( A )
820    RETURN
      END
```




Figure 18. Flowchart of Subroutine CALC

```

SUBROUTINE CALC

C      THIS SUBROUTIN READS IN THE DESIRED DATA FILES
C      AND CALLS CRUNCH WHICH PERFORMS THE
C      ACTUAL CALCULATIONS

C      * * * COMMON VARIABLES * * *
      LOGICAL*1      HAZ(3),HAZA(5)
      INTEGER        I2(256)
      INTEGER*4      I4(12)
      REAL*4         R4(64)
      INTEGER        IDATA(1000)

      COMMON /HAZEL/HAZ,HAZA
      COMMON /VARBLE/I2,I4,R4
      COMMON /TEMPS/IDATA
C      * * * END COMMON BLOCK * * *

C      * * * BEGIN COMMON DESIGNATIONS * * *
      INTEGER IACTVS(6),IACTVC(6),IPROBE(6),IHEAT(6),IPRBS9(6)
      INTEGER IFILE(6,3),ITEXT(6,20),ICDATE(6,3)

      EQUIVALENCE      (ICDATE(1,1),I2(150)),      (IACTVC(1),I2(168))
      EQUIVALENCE      (IFILE(1,1),I2(174)),      (IAUTOC,I2(192))
C      * * * END COMMON DESIGNATIONS * * *

C      * * * BEGIN LOCAL VARIABLES * * *
      INTEGER IFILEA(7)
C      * * * END LOCAL VARIABLES * * *

C      * * * BEGIN CODE SEGMENT * * *
      IFILEA(1)='FD'
      IFILEA(2)='1'
      IFILEA(6)='D'
      IFILEA(7)='AT'

      HAZ(1)=26
      CALL PRINT(HAZ)

      DO 100 I=1,6
      IF(IACTVC(I).EQ.0) GOTO 100
      TYPE 5,I
5      FORMAT('CRUNCH CHANNEL NUMBER ',I)

      DO 10 J=1,3
      IFILEA(J+2)=IFILE(I,J)
10     CONTINUE

      CALL ASSIGN(3,IFILEA,14)
      READ(3,ERR=60) IDATA
60     CALL CLOSE(3)

      CALL CRUNCH

      TYPE 70
70     FORMAT('BACK FROM CRUNCH')

100    CONTINUE
      ACCEPT 30
      FORMAT( A )
      RETURN
      END

```





Figure 19. Flowchart of Subroutine AUTO

```

SUBROUTINE AUTO

C      THIS SUBROUTINE CALL SAMPLE AND CALCULATE,
C      INCREMENTS FILENAMES AND ADJUSTS
C      REPETITION PARAMETERS

C      * * * COMMON VARIABLES * * *
      LOGICAL*1      HAZ(3),HAZA(5)
      INTEGER         I3(256)
      INTEGER*4       I4(12)
      REAL*4          R4(64)
      INTEGER         IDATA(1800)

      COMMON /HAZEL/HAZ,HAZA
      COMMON /VARBLE/I2,I4,R4
      COMMON /TEMPS/IDATA
C      * * * END COMMON BLOCK * * *

C      * * * BEGIN COMMON DESIGNATIONS * * *
      INTEGER IACTVS(6),IACTVC(6),IPROBE(6),IHEAT(6),IPPRBB(6)
      INTEGER IFILE(6,3),ITEXT(6,20),ICDATE(6,3)

      EQUIVALENCE      (ICOM,I2(1)),          (IACTVS(1),I2(2))
      EQUIVALENCE      (ICDATE(1,1),I2(150)), (IACTVC(1),I2(168))
      EQUIVALENCE      (IFILE(1,1),I2(174)),   (IAUTOC,I2(192))
      EQUIVALENCE      (IEXREP,I2(193))
      EQUIVALENCE      (RLALIM,R4(44)),        (AUTOTN,R4(45))
      EQUIVALENCE      (AUTOTU,R4(46)),        (EXINT,R4(47))
C      * * * END COMON DESIGNATIONS * * *

C      * * * BEGIN LOCAL VARIABLES * * *
      INTEGER*4 KTIME
      INTEGER    KTIM(2)

      EQUIVALENCE (KTIME,KTIM(1))
C      * * * END LOCAL VARIABLES * * *

C      * * * BEGIN CODE SEGMENT * * *
      CALL GTIM(KTIME)
      AUTOTU=65536.0*KTIM(1)+KTIM(2)
      AUTOTN=EXINT

C      * * * PERFORM SAMPLE (IF NEEDED) * * *
      J=0
      DO 100 I=1,6
      IF(IACTVS(I).NE.0) J=J+1
100    CONTINUE
      IF(J.NE.0) CALL SAMP

C      * * * PERFORM CALCULATE (IF NEEDED) * * *
      J=0
      DO 110 I=1,6
      IF (IACTVC(I).NE.0) J=J+1
110    CONTINUE
      IF (J.NE.0) CALL CALC

C      * * * INCREMENT FILENAMES * * *
      DO 120 I=1,6
      IF(IACTVS(I).EQ.0.AND. IACTVC(I).EQ.0) GOTO 180
      DECODE(2,129,IFILE(I,3))IFILEA
120    FORMAT(I2)
      IFILEA=IFILEA+1
      TYPE 125,I,IFILE(I,3),IFILEA
125    FORMAT(1,' ',A2,' ',I3)
      IF(IFILEA.LT.100) GOTO 160

```

```

      IFILEA=0
C   * * *   ADD 120 * * *
      IFILES=IFILE(I,2)+*100
      TYPE 140, IFILE(I,2), IFILES, IFILES
140   FORMAT(A2, ' ', A2, ' ', 07)
C   * * *   CHECK GE TO XXX900 * * *
C   * * *   LEAVE AT XXX899 * * *
      IF(IFILES.GE."34400") GOTO 180
      IFILE(I,2)=IFILES
150   ENCODE(2,170, IFILE(I,3)) IFILEA
170   FORMAT(I2)
      IF(IFILEA.LT.10) IFILE(I,3)=IFILE(I,3).OR."60
      IF(IFILEA.EQ.0) IFILE(I,3)="30060
180   CONTINUE

C   * * *   UPDATE AUTO PARAMETERS * * *
      IEXREP=IEXREP-1
      IF(IEXREP.NE.0) GOTO 200
      ICOM='B'
      AUTOTN=0.0
      IEXREP=1

200   RETURN
      END

```




Figure 20. Flowchart of Subroutine DATA

SUBROUTINE DATA

```

C      THIS SUBROUTINE ALLOWS THE USER TO EITHER LIST
C      AN ENTIRE DATA FILE OR PLOT A FILE ON THE
C      GRAPHICS TERMINAL

C      * * * COMMON VARIABLES * * *
      LOGICAL*1      HAZ(3),HAZA(5)
      INTEGER        I2(256)
      INTEGER*4      I4(12)
      REAL*4         R4(64)
      INTEGER        IDATA(1800)

      COMMON /HAZEL/HAZ,HAZA
      COMMON /VARBLE/I2,I4,R4
      COMMON /TEMPS/IDATA
C      * * * END COMMON BLOCK * * *

C      * * * BEGIN LOCAL DESIGNATIONS * * *
      INTEGER*4 ISTIM
      INTEGER KDATE(3),IFILET(18),ITEXT(20)
      INTEGER ICOATE(3)
      REAL*4 PPR80,PPBETA,PPA,PPS,PPV88
      INTEGER*4 IPPT88
      REAL*4 DATA(1800)

      INTEGER IFILEA(7)
      COMMON /GRID/IX,IY
C      * * * END LOCAL DECLARATIONS * * *

      IFILEA(1)='FD'
      IFILEA(2)='1.'
      IFILEA(6)='D'
      IFILEA(7)='AT'
10      JTEMP=ITTOUR(26)
      TYPE 30
30      FORMAT( )
      TYPE 40
40      FORMAT( 'ENTER FILENAME YOU WISH TO OPEN  ',*)
      ACCEPT 50,IFILEA(3),IFILEA(4),IFILEA(5)
50      FORMAT( A2,A2,A2 )
      CALL ASSIGN(3,IFILEA,14)
      READ (3,ERR=60) J,ISTIM,KDATE,I,IFILET,ITEXT,IProbe,
1      PPR80,PPBETA,PPA,PPS,ICOATE,IPPR88,IPPT88,PPV88,
2      IHEAT,IPSR,IPSP,ISR,ISP,IDATA
60      CALL CLOSE(3)
      ICHAN=I

65      TYPE 70
70      FORMAT( 'ENTER "G" FOR GRAPHICS OR "L" FOR LISTING')
      ACCEPT 80,I
80      FORMAT(A)
      IF(I.EQ.'G') GOTO 500

C      * * * LISTING PACKAGE * * *
      I=ICHAN
      TYPE 100,J
100     FORMAT('DATA VERSION - - ',I)
      CALL CVTTIM(ISTIM,KHRS,KMIN,KSEC,KTICKS)
      TYPE 110,KHRS,KMIN,KSEC,KDATE
110     FORMAT('SAMPLE TIME - - ',I2,'.',I2,'.',I2,
1      ' DATE - - ',I2,'.',I2,'.',I2)
      TYPE 120,I
120     FORMAT('CHANNEL DATA RECIEVED FROM - - ',I)

```

```

TYPE 130
FORMAT('ACTIVE FILES AT SAMPLE TIME')
TYPE 140,IFILET
140 FORMAT('1- ',3A2,' 2- ',3A2,' 3- ',3A2,' 4- ',3A2,' 5- ',
1      3A2,' 6- ',3A2)
TYPE 150,ITEXT
150 FORMAT('TEXT - - ',20A2)
TYPE 155,IPROBE
155 FORMAT('PROBE DESIGNATION - - ',I3)
TYPE 160,PPR80,PPBETA,PPA,PP8
160 FORMAT('R80 - - ',F7.2,' BETA - - ',F8.3,
1      ' A - - ',F8.4,' B - - ',F8.4)
TYPE 165,ICDATE
165 FORMAT('CALIBRATED ON - ',I2,'/',I2,'/',I2)
CALL CVTTIM(IPPT88,KHRS,KMIN,KSEC,KTICKS)
TYPE 170,IPPR88,PPV88,KHRS,KMIN,KSEC
170 FORMAT('BRIDGE BALANCE : OHMS - - ',I6,' VOLTS - - ',F5.2,
1      ' TIME - - ',I2,' ',I2,' ',I2)
TYPE 180,IHEAT
180 FORMAT('HEAT PULSE DURATION - - ',I4)
TYPE 190,IPSR,IPSP
190 FORMAT('PRESAMPLE READINGS - - ',I4,
1      ' PRESAMPLE PERIOD - - ',I3)
TYPE 200,ISR,ISP
200 FORMAT(' SAMPLE READINGS - - ',I4,
1      ' SAMPLE PERIOD - - ',I3)
PPR80=5.25/'3777
DO 205 K=1,IPSR+ISR
IF(IDATA(K).GT.'3777') IDATA(K)=IDATA(K)-'10000
DATA(K)=IDATA(K)*PPR80
205 CONTINUE
TYPE 210
210 FORMAT('PRESAMPLE DATA : (MILLIVOLTS)')
DO 225 I=1,IPSR,10
J=I+9
IF(IPSR-I.LT.9) J=IPSR
TYPE 220,(DATA(K),K=I,J)
FORMAT(10F7.4)
225 CONTINUE
TYPE 230
230 FORMAT('SAMPLE DATA : (MILLIVOLTS)')
DO 245 I=IPSR+1,IPSR+ISR,10
J=I+9
IF(IPSR+ISR-I.LT.9) J=IPSR+ISR
TYPE 240,(DATA(K),K=I,J)
FORMAT(10F7.4)
245 CONTINUE
GOTO 290

C * * * GRAPHICS PACKAGE * * *
500 JTEMP=ITTOUR(25)

C * * * SET GRAPH SIZE
JTEMP=ITTOUR(13)
IXMIN=35
IXMAX=1010
IYMIN=40
IYMAX=740
TYPE 510
510 FORMAT('DO YOU WANT DIFFERENT SIZE GRAPH?')
ACCEPT 520,I
520 FORMAT(A)
IF(I NE. 'Y') GOTO 545
TYPE 530
530 FORMAT('CONSIDER THE SCREEN TO BE 100 BY 100')

```

```

332      TYPE 935
335      FORMAT('ENTER LEFT HORIZONTAL MARGIN FROM 0 TO 99')
      ACCEPT 940,I
340      FORMAT(I2)
      TYPE 945,I
345      FORMAT('ENTER RIGHT HORIZONTAL MARGIN FROM ',I2,' TO 99')
      ACCEPT 950,J
350      FORMAT(I2)
      IF(J.GT.I)GOTO 960
      TYPE 955
355      FORMAT('INCORRECT VALUES - - TRY AGAIN')
      GOTO 932
360      I=<(IXMAX-IXMIN)*I/99+IXMIN
      IXMAX=INT<<(IXMAX-IXMIN)/99.0*J>+IXMIN
      IXMIN=I

362      TYPE 965
365      FORMAT('ENTER LOWER VERTICAL LIMIT FROM 0 TO 99')
      ACCEPT 970,I
370      FORMAT(I2)
      TYPE 975,I
375      FORMAT('ENTER UPPER VERTICAL LIMIT FROM ',I2,' TO 99')
      ACCEPT 980,J
380      FORMAT(I2)
      IF(J.GT.I)GOTO 990
      TYPE 985
385      FORMAT('INCORRECT VERTICAL LIMITS - - TRY AGAIN')
      GOTO 962
390      I=<(IYMAX-IYMIN)*I/99+IYMIN
      IYMAX=<(IYMAX-IYMIN)*J/99+IYMIN
      IYMIN=I

345      JTEMP=ITTOUR(26)

C   * * *   SET TIME SCALE LIMITS   * * *
      ITNEG=IPSR*IPSP+IHEAT
      ITPOS=ISR*ISP
      ITMINS=-ITNEG/60
      ITMAXS=ITPOS/60
      ISTART=ITNEG+ITMINS*60
      ITOTAL=ITNEG+ITPOS

      TYPE 1010
1010     FORMAT('DO YOU WANT TO EXPAND TIME SCALE?')
      ACCEPT 1015,I
1015     FORMAT(A)
      IF(I.NE.'Y') GOTO 1100
      TYPE 1025,ITMINS,ITMAXS
1020     FORMAT('ENTER LOWER TIME BOUNDARY - - FROM ',I3,' TO ',I3)
1025     ACCEPT 1027,I
1027     FORMAT(I3)
      IF(I.GE.ITMAXS.OR.I.LT.ITMINS) GOTO 1020
1029     TYPE 1030,I,ITMAXS
1030     FORMAT('ENTER UPPER TIME BOUNDARY - - FROM ',I3,' TO ',I3)
      ACCEPT 1032,J
1032     FORMAT(I3)
      IF(J.GT.ITMAXS.OR.J.LE.I) GOTO 1029
      ITMAXS=J
      ITMINS=I
      ISTART=0
      ITOTAL=<(ITMAXS-ITMINS)*60
1100     JTEMP=ITTOUR(25)

C   * * *   LABEL CURVE   * * *

```

```

JTEMP=ITTOUR(29)
IX=0
IY=779
CALL GPOINT
JTEMP=ITTOUR(13)
TYPE 810, ITEXT
310 FORMAT(30X, 20A2)
TYPE 820, IFILEA(3), IFILEA(4), IFILEA(5), ICHAN
820 FORMAT(30X, 'FILENAME: ', 3A2, ' CHANNEL: ', I)
TYPE 830, IPROBE, IHEAT/60, IHEAT-IHEAT/60*60
830 FORMAT(30X, 'PROBE: ', I2, ' HEAT PULSE: ', I2, ' ', I2)
CALL CVTTIM(ISTIM, KHRS, KMIN, KSEC, KTICKS)
TYPE 840, KHRS, KMIN, KSEC, KOATE
840 FORMAT(30X, 'SAMPLE TIME- ', I2, ' ', I2, ' ', I2, ' DATE- ',
1 I2, ' ', I2, ' ', I2, ' ', I2)
JTEMP=ITTOUR(28)

```

```

C * * * DRAW AXIS * * *
JTEMP=ITTOUR(29)
IX=IXMIN
IY=IYMAX
CALL GPOINT
IX=IXMIN
IY=IYMIN
CALL GPOINT
IX=IXMAX
IY=IYMIN
CALL GPOINT
JTEMP=ITTOUR(28)

```

```

C * * * LABEL AXIS * * *
B=(IYMAX-IYMIN)/10.5
YMID=(IYMAX+IYMIN)/2.0+IYMIN
JTEMP=ITTOUR(31)
TYPE 548
548 FORMAT( )
DO 580 I=-5.5
IY=INT(YMID+I*8) + 11
IX=IXMIN-35
JTEMP=ITTOUR(29)
CALL GPOINT
JTEMP=ITTOUR(31)
TYPE 550, I
550 FORMAT( I2)
JTEMP=ITTOUR(29)
IY=IY-11
IX=IXMIN-5
CALL GPOINT
IX=IXMIN
CALL GPOINT
JTEMP=ITTOUR(28)
558 CONTINUE

D=FLOAT(IXMAX-IXMIN)/ITOTAL
DO 630 I=ITMINS, ITMAXS
IX=INT((ISTART+(I-ITMINS)*60)*D)+IXMIN
IY=IYMIN
JTEMP=ITTOUR(29)
CALL GPOINT
IY=IY-5
CALL GPOINT
JTEMP=ITTOUR(28)
JTEMP=ITTOUR(29)

```



```

IX=IX-15
CALL GPOINT
JTEMP=ITTOUR(31)
TYPE 640,I
540  FORMAT(I2)
JTEMP=ITTOUR(28)
650  CONTINUE

C * * * PLOT CURVE * * *
ICOUNT=0
ITLOW=ITMINS*60+ITNEG-ISTART
ITHIGH=ITMAXS*60+ITNEG-ISTART
D=FLOAT(IXMAX-IXMIN)/(ITHIGH-ITLOW)
B=FLOAT(IYMAX-IYMIN)/"10000
DO 720 I=1,IPSR+ISR
IT=IPSP*I
IF(I.GT.IPSR) IT=ITNEG+ISP*(I-IPSR)
IF(IT.LT.ITLOW.OR.IT.GT.ITHIGH) GOTO 720
IF(IT.LT.ITNEG.AND.IT.GT.IPSP*IPSR) GOTO 720
IX=IXMIN+INT((IT-ITLOW)*D)
ICOUNT=ICOUNT+1
IF(ICOUNT.NE.15)GOTO 705
ICOUNT=0
JTEMP=ITTOUR(31)
TYPE 703
703  FORMAT( )
JTEMP=ITTOUR(28)
705  IY=IDATA(I)
IF(IY.GT."3777") IY=IY-"10000
IY=INT(YMID+8*IY)
CALL GPOINT
720  CONTINUE

JTEMP=ITTOUR(31)
JTEMP=ITTOUR(24)
HAZA(3)=37
HAZA(4)=72
CALL PRINT(HAZA)

C * * * RAP THINGS UP * * *
290  TYPE 295
295  FORMAT('ENTER "E" TO EXIT OR "D" FOR MORE DATA ',S)
ACCEPT 300,I
300  FORMAT(A)
JTEMP=ITTOUR(28)
JTEMP=ITTOUR(31)
JTEMP=ITTOUR(25)
JTEMP=ITTOUR(24)
JTEMP=ITTOUR(26)
IF(I.EQ.'E') GOTO 310
TYPE 303
303  FORMAT('DO YOU WANT TO REMAIN WITH SAME DATA ?')
ACCEPT 307,I
307  FORMAT(A)
IF(I.EQ.'Y') GOTO 65
GOTO 10
310  RETURN
END

```




Figure 21. Flowchart of PROBEC.FOR

```

1      PROSEC FOR
2
3      THIS PROGRAM MODIFIES THE FILE PROBE.DAT
4      AND USES SUBROUTINE GET TO ACCEPT NEW CALIBRATIONS
5
6      DIMENSION IDATA(1500)
7      COMMON IDESIG,IMONTH,IDAY,IYEAR,ATEM,BTEM,RSOTEM,BETAT
8      DIMENSION JDATA(12)
9      EQUIVALENCE (ATEM,JDATA(5)),(BTEM,JDATA(7))
10     EQUIVALENCE (RSOTEM,JDATA(9)),(BETAT,JDATA(11))
11     EQUIVALENCE (IDESIG,JDATA(1)),(IMONTH,JDATA(2))
12     EQUIVALENCE (IDAY,JDATA(3)),(IYEAR,JDATA(4))
13     CALL ASSIGN(3,'F00:PROBE.DAT',13)
14     READ (3) JNUM
15     INUM=JNUM*12
16     READ (3) (IDATA(K),K=1,INUM)
17     CALL CLOSE(3)
18
19     TYPE 6
20     FORMAT( )
21     TYPE 2
22     FORMAT( 'ENTER "1" TO ADD OR MODIFY A PROBE ' )
23     TYPE 3
24     FORMAT( 'ENTER "2" TO LIST ALL PROBES' )
25     TYPE 4
26     FORMAT( 'ENTER "3" TO EXIT' )
27     ACCEPT 8,ICOM
28     FORMAT( I )
29     IF (ICOM.EQ.2) GOTO 200
30     IF (ICOM.EQ.3) GOTO 300
31     IF (ICOM.NE.1) GOTO 5
32     TYPE 6
33     TYPE 25
34
35 C * * * ADD OR MODIFY A PROBE * * *
36 25     FORMAT( 'ENTER PROBE NUMBER' )
37     ACCEPT 30,IDESIG
38     FORMAT( I3 )
39     TYPE 6
40     INDEX=1
41     IF (IDESIG.EQ.IDATA(INDEX)) GOTO 150
42     INDEX=INDEX+12
43     IF (INDEX.LT.INUM) GOTO 40
44     TYPE 100
45     FORMAT( 'NEW PROBE NUMBER' )
46     TYPE 110
47     FORMAT( 'ENTER "Y" TO START A NEW PROBE WITH THAT' )
48     1     ' DESIGNATION' )
49     ACCEPT 112,ICOM
50     FORMAT( A )
51     TYPE 6
52     IF (ICOM.NE.'Y') GOTO 5
53     CALL IDATE,IMONTH,IDAY,IYEAR)
54     ATEM=0
55     BTEM=0
56     RSOTEM=0
57     BETAT=0
58     JNUM=JNUM+1
59     INDEX=INUM+1
60     INUM=INUM+12
61     CALL GET
62     DO 122 K=1,12
63     IDATA(INDEX)=JDATA(K)
64     INDEX=INDEX+1
65 122     CONTINUE

```

```

      GOTO 5
150  TYPE 155
155  FORMAT( 'PROBE CURRENTLY HAS THESE VALUES' )
      DO 150 K=1,12
      JDATA(K)=IDATA(INDEX)
      INDEX=INDEX+1
160  CONTINUE
      INDEX=INDEX-12
      CALL GET
      GOTO 129
200  TYPE 6
      TYPE 202

C * * * LIST ALL PROBES IN FILE * * *
332  FORMAT( 'THIS IS A COMPLETE LISTING OF PROBES AND DATES' )
      DO 210 K=1,JNUM
      INDEX=K*12-11
      DO 205 J=1,4
      JDATA(J)=IDATA(INDEX)
      INDEX=INDEX+1
235  CONTINUE
      TYPE 207, IDESIG, IMONTH, IDAY, IYEAR
207  FORMAT( I3, 5X, I2, '/', I2, '/', I2 )
210  CONTINUE
      GOTO 5

C * * * CHECK IF REPLACEMENT DESIRED * * *
330  TYPE 310
310  FORMAT( 'LAST CHANCE TO REMAIN WITH EXISTING PROBE FILE' )
      TYPE 315
315  FORMAT( 'ENTER "1" TO REPLACE OLD PROBE FILE' )
      ACCEPT 320, ICOM
320  FORMAT( I )
      IF (ICOM.NE.1) GOTO 400
      CALL ASSIGN(3, 'F00', 'PROBE.DAT', 13)
      WRITE(3) JNUM
      WRITE(3) IDATA
      CALL CLOSE(3)
      TYPE 325
325  FORMAT( 'NEW FILE CREATED - - OLD FILE DESTROYED' )
      GOTO 490
400  TYPE 410
410  FORMAT( 'NO CHANGES MADE TO EXISTING PROBE FILE' )
490  FORMAT( 'STILL UNDER CONSTRUCTION' )
490  STOP ' THAT IS ALL FOR THIS MESS '
      END

C          CREATE PROBE.DAT
C
C  THIS PROGRAM CREATES AN EMPTY FILE
C  TO BE USED BY PROBE.C.FOR
C
      DIMENSION JDATA(12)
      CALL ASSIGN(3, 'F00', 'PROBE.DAT', 13)
      WRITE(3) 1
      DO 5 K=1,12
      JDATA(K)=0
5      CONTINUE
      WRITE(3) JDATA
      CALL CLOSE(3)
      STOP 'FILE CREATED'
      END

```




Figure 22. Flowchart of Subroutine GET

SUBROUTINE GET

```

0      THIS SUBROUTINE ALLOWS THE USER TO CHANGE THE
0      CALIBRATION DATA FOR A GIVEN PROBE

      COMMON IDESIG,IMONTH,IDAY,IYEAR,ATEM,BTEM,RB0TEM,BETAT
7      TYPE 5,IDESIG
8      FORMAT( 'PROBE NUMBER = ',I3,' COMMAND' )
      TYPE 10,IMONTH
10     FORMAT( 'MONTH      = ',I3,' (1)' )
      TYPE 15,IDAY
15     FORMAT( 'DAY        = ',I3,' (2)' )
      TYPE 20,IYEAR
20     FORMAT( 'YEAR       = ',I3,' (3)' )
      TYPE 25,ATEM
25     FORMAT( 'A          = ',F8.4,' (4)' )
      TYPE 30,BTEM
30     FORMAT( 'B          = ',F8.4,' (5)' )
      TYPE 35,RB0TEM
35     FORMAT( 'RB0       = ',F8.2,' (6)' )
      TYPE 40,BETAT
40     FORMAT( 'BETA      = ',F8.3,' (7)' )
      TYPE 45
45     FORMAT( 'TO EXIT TYPE      (8)' )
      ACCEPT 100,ICOM
      FORMAT( I )
100    IF (ICOM.NE.1) GOTO 110
      TYPE 105
105    FORMAT( 'ENTER MONTH' )
      ACCEPT 106,IMONTH
106    FORMAT( I2 )
110    IF (ICOM.NE.2) GOTO 120
      TYPE 115
115    FORMAT( 'ENTER DAY' )
      ACCEPT 106,IDAY
120    IF (ICOM.NE.3) GOTO 130
      TYPE 125
125    FORMAT( 'ENTER YEAR' )
      ACCEPT 106,IYEAR
130    IF (ICOM.NE.4) GOTO 140
      TYPE 135
135    FORMAT( 'ENTER A' )
      ACCEPT 136,ATEM
136    FORMAT( F8.4 )
140    IF (ICOM.NE.5) GOTO 150
      TYPE 145
145    FORMAT( 'ENTER B' )
      ACCEPT 136,BTEM
150    IF (ICOM.NE.6) GOTO 160
      TYPE 155
155    FORMAT( 'ENTER RB0' )
      ACCEPT 136,RB0TEM
160    IF (ICOM.NE.7) GOTO 170
      TYPE 165
165    FORMAT( 'ENTER BETA' )
      ACCEPT 136,BETAT
170    IF (ICOM.NE.8) GOTO 3
      RETURN
      END

```

REFERENCES

1. Chen, M. M. and K. R. Holmes. "The thermal pulse-decay method for simultaneous measurement of thermal conductivity and local blood perfusion rate of living tissues." 1980 Advances in Bioengineering, Am. Soc. Mech. Engr., 1980, pp. 113-115.
2. Aukland, K. "Methods for measuring renal blood flow: total flow and regional distribution." Ann. Rev. Physiol. 42:543-555, 1980.
3. Holmes, K. R. and M. M. Chen. "In vivo tissue thermal conductivity and local blood perfusion measured with a heat pulse-decay method." Physiologist 23:183, 1980.
4. Holmes, K. R. and M. M. Chen. Local thermal conductivity of Para-7 fibrosarcoma in hamster. 1979 Advances in Bioengineering, Am Soc. Mech. Engr., 1979, 147-149.
5. Digital Equipment Corp. Microcomputer Processors (1978).
6. Digital Equipment Corp. Memories and Peripherals (1979).
7. Adac Corporation. Model 1030 Analog to Digital Converter Manual (1981).
8. Lear Sigler, Inc. ADM 5 Dumb Terminal Video Display Unit Users Reference Manual (1981).
9. Digital Engineering, Inc. User's Manual RG-512 Retro-graphics Card for the ADM-3A Computer Terminal (1981).
10. Digital Engineering, Inc. GP-100 Graphx Printer User's Manual (1981).
11. Digital Equipment Corp. RT-11 System Reference Manual (1979).